# Symphony®

**EN** Breast pump –  
for healthcare professionals  
Instructions for use

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# 1. Important safety information

4

Any serious incident that has occurred in relation to the Symphony breast pump should be reported to Medela AG and the relevant Competent Authority.

**Please read all instructions before using this product.  
Keep these instructions for future reference.**

Warnings identify all instructions that are important for safety. Failure to observe these instructions can lead to personal injury or damage to the product.

The following symbols and signal words show the significance of the warnings:

 **WARNING**

Can lead to serious injury or death.

 **CAUTION**

Can lead to minor injury.

**NOTICE**

Can lead to material damage  
(not related to personal injury).

 **Information**

Useful or important information  
that is not related to safety.

 **WARNING**

To reduce the risk of electric shock:

- Plug the power cord into the breast pump first and then into the wall socket.
- Keep the device dry. Never immerse in water or other liquids.
- Do not spray or pour liquid directly onto the pump unit.
- Do not place or store the device where it can fall or be pulled into a sink or tub.
- Do not use the Symphony breast pump while bathing or showering.
- Do not reach for any electrical device if it has fallen into water. Unplug immediately.

 **WARNING**

To avoid fire, electric shock or serious burns:

- The Symphony breast pump is not heat-resistant: keep away from radiators and open flames.
- Do not expose the pump unit to direct sunlight.
- Do not operate the Symphony breast pump where aerosol (spray) products are being used or where oxygen is being administered.
- Do not use the device in the presence of a flammable anaesthetic mixture with air, or with oxygen or nitrous oxide.
- Do not thaw frozen breast milk, or do not warm breast milk in a microwave or in a pan of boiling water. Microwaving can cause severe burns to the baby's mouth from hot spots that develop in the milk during microwaving. Microwaving can also change the composition of breast milk.
- Only use the power cord that comes with the Symphony breast pump.
- Never operate an electrical device if it has a damaged power cord, if it is not working properly, if it has been dropped or damaged, or dropped into water.
- Always inspect power cord wires prior to use for damage or exposed wire. If any damage is found, immediately stop using the power cord. Contact the local Medela representative.
- The Symphony breast pump should never be left unattended when plugged into a power source.
- Separation from the mains is only assured by disconnecting the power cord from the wall socket.
- Position the pump unit such that it is easy to disconnect the power cord from the wall socket.
- Keep the power cord away from heated surfaces.

## Important safety information (cont.)

6

### **WARNING**

To avoid health risk and reduce risk of injury:

- Do not continue pumping for more than two consecutive pumping sessions if no results are achieved.
- If problems or pain occur, consult a lactation consultant or breastfeeding specialist.
- Never use while pregnant, unless prescribed/advised by a physician, as pumping can induce labour.
- For mothers who are infected with Hepatitis B, Hepatitis C or Human Immuno-deficiency Virus (HIV): pumping breast milk will not reduce or remove the risk of transmitting the virus to the baby through breast milk.
- Inspect all components of the Symphony breast pump as well as all parts of the pump set(s) before each use.
- Never use a device, components or parts that are damaged.
- Never use a device, components or parts if dirt, mould or other contamination is visible.
- Replace mouldy, damaged or worn parts.
- There are no user-serviceable parts inside the Symphony pump unit. Repairs must only be performed by an authorised service agency. Do not attempt to repair the pump or parts by yourself.
- Only use the Symphony breast pump for its intended use as described in these instructions for use.
- Only use Medela original accessories.
- Do not modify any parts of the pump or the pump set.
- Do not use the Symphony breast pump while operating a motor vehicle.
- Do not drive whilst hands-free pumping.
- Never drop or insert any object into any opening or tubing.
- Do not use outdoors.
- Supervision is necessary when the Symphony breast pump is used in the vicinity of children.
- Portable and mobile radio frequency communications equipment can affect the Symphony breast pump.
- The Symphony breast pump needs special precautions regarding EMC (Electro Magnetic Compatibility). Install and put into service according to the EMC information provided in these instructions for use.

**⚠ CAUTION**

To avoid risk of contamination:

- Wash hands thoroughly with soap and water before touching the breast pump, the pump set(s) and breasts. Avoid touching the inside of containers or lids.
- Do not use the Symphony breast pump while sleeping or overly drowsy.
- Only use drinking-quality tap or bottled water for cleaning.
- Do not store wet or damp parts as mould may develop.

To avoid health risk and reduce risk of injury:

- Only pump with the breast pump in an upright position.
- Pressing shields too firmly onto the breasts may affect the milk flow.
- Do not express with a vacuum that is too high and uncomfortable (painful). Pain, along with potential breast and nipple trauma, may decrease the milk output.
- Do not attempt to remove the breast shield from the breast while pumping. Turn the breast pump off and break the seal between the breast and the breast shield with a finger, then remove the breast shield from the breast.
- If pumping is uncomfortable or causing pain, turn the unit off, break the seal between breast and breast shield with a finger and remove the breast shield from the breast.
- Do not shorten the tubing between breast shield and pump.
- Do not use abrasive washing-up liquids/detergents when cleaning the breast pump, pump components or pump set parts.
- Never put the pump unit in water or in a steriliser, as this can cause permanent damage to the breast pump.
- Do not use Medela Symphony pump sets with non-compatible breast pumps.
- Do not use the Medela Symphony pump sets for any other purpose than expressing breast milk.

**NOTICE**

- Plastic bottles and component parts become brittle when frozen and may break when dropped.
- Also, bottles and component parts may become damaged if mishandled, e.g. dropped, over-tightened, or knocked over.
- Take appropriate care in handling bottles and components.
- Do not use the breast milk if bottles or components show signs of damage.

<b>1. Important safety information</b>	<b>4</b>
<b>2. Intended use/indications</b>	<b>10</b>
2.1 Intended purpose	10
2.2 Indications for use	10
2.3 Contraindications	10
<b>3. Product description</b>	<b>10</b>
3.1 Symphony hospital breast pump and accessories	10
3.2 Symphony hospital breast pump – product parts	12
3.3 Reusable pump set – product parts (applied parts)	13
3.4 Disposable pump set – product parts (applied parts)	13
<b>4. Cleaning</b>	<b>14</b>
4.1 Breast pump – cleaning and disinfection	14
4.2 Cleaning pump sets	14
<b>5. Preparing for pumping</b>	<b>15</b>
5.1 Installation of the program card	15
5.2 Installation of the card and cord protector (optional)	15
5.3 Mains operation	16
5.4 Assembling the Symphony pump set	16
5.5 Connecting the pump set to the pump housing	17
<b>6. Expression</b>	<b>18</b>
6.1 Selecting the correct Medela breast shield size	18
6.2 Selecting the right program	20
6.3 Operation: INITIATE program	21
6.4 Operation: MAINTAIN program	22
6.5 Set maximum comfort vacuum	24
6.6 Expression by single pumping	24
6.7 Expression by double pumping	25
<b>7. Storing and thawing breast milk</b>	<b>26</b>
7.1 Storage and thawing guidelines	26

<b>8. Troubleshooting and error codes</b>	<b>27</b>
8.1 Troubleshooting	27
8.2 Error codes	28
<b>9. Warranty and maintenance</b>	<b>29</b>
<b>10. Disposal</b>	<b>30</b>
<b>11. Accessories list</b>	<b>31</b>
<b>12. Meaning of symbols</b>	<b>34</b>
<b>13. Electromagnetic compatibility (EMC)</b>	<b>36</b>
<b>14. Technical specifications</b>	<b>37</b>
<b>15. Appendix</b>	<b>208</b>
15.1 EMC/technical description (English only)	208

## 2. Intended use/indications

10

### 2.1 Intended purpose

The Symphony breast pump is a powered breast pump to be used by lactating women to express and collect milk from their breasts.

### 2.2 Indications for use

The Symphony breast pump is indicated to express and collect milk from the breasts of a lactating woman. It is also indicated to initiate, build and maintain a mother's milk supply to support the provision of human milk for infants who are not able to directly or effectively breastfeed (e.g. latch-on problems such as flat or inverted nipples, cleft palate and prematurity). The breast pump removes milk from affected breast/s to alleviate the symptoms of engorgement and milk stasis as well as support the healing process of sore or cracked nipples and mastitis.

### 2.3 Contraindications

There are no known contraindications for the Symphony breast pump.

## 3. Product description

### 3.1 Symphony hospital breast pump and accessories

Medela's Symphony hospital-grade breast pump is a multi-user breast pump with research-based pumping programs which were developed specifically to support mothers of preterm and term infants to initiate, build and maintain an adequate milk production.

For adequate initiation and building of milk supply, the breasts should be drained

- early: within 1-3 hours after birth
- frequently: every 2 -3 hours
- effectively: baby with a healthy sucking behaviour or double pumping.

It is recommended to record the date, time and volumes pumped.

The Symphony breast pump is provided with the Symphony PLUS program card, which contains two pump suction pattern programs.

### **The INITIATE program with Medela's initiation technology**

The INITIATE program mimics the irregular and more rapid sucking and pausing pattern of a term born infant in the first few days of lactation. This program is intended for use in the first few days after birth before secretory activation to support pump-dependent mothers to initiate milk production.

### **The MAINTAIN program with Medela's 2-Phase Expression technology**

The MAINTAIN program mimics the term born infant's sucking pattern during established lactation. This program is designed to optimise milk output after secretory activation to support all mothers to build and maintain lactation.

In some cases, the Symphony breast pump is provided with the Standard program card. This program card contains the MAINTAIN program. While the Standard card can be used to initiate lactation, Medela recommends using the INITIATE program of the Symphony PLUS program card.

## **Required accessories**

Only operate the Symphony breast pump with Medela pump sets and milk bottles. Pump sets and bottles are available as reusable or as disposable products.

### **Reusable products:**

Symphony Hospital Multibox (25 reusable pump sets in sizes S or M).

### **Disposable products:**

Medela's disposable products are intended for single-use or one-day use in hospitals, without a need for cleaning prior to first use:

Medela's Ready-to-Use products are hygienically safe to use without prior cleaning, for all hospitals where the use of non-sterile products is allowed.

Medela's sterile products are EO sterile and fulfil the EN/ISO standards for medical devices ISO 11607-1, ISO 11607-2, ISO 11135.

# Product description (cont.)


## Available disposable products:

- Sterile single-use Symphony Pump Set (size S, M, L or XL).
- Ready-to-Use one-day Symphony Pump Set (size S, M, L or XL).
- Colostrum container (35 ml).
- Sterile single-use bottle (available sizes: 80, 150 and 250 ml).
- Ready-to-Use single-use bottle (available sizes: 80, 150 and 250 ml).

## Optional accessories and other products to be used in combination

Mobile stand for the breast pump, hard case for the pump unit, card and cord protector.


### 3.2 Symphony hospital breast pump – product parts

Accessories available for order can be found in **chapter 11**. Articles may not be available in all countries. For further information about Medela products, visit [www.medela.com](http://www.medela.com) and look for location finder.


### 3.3 Reusable pump set – product parts (applied parts)

Symphony reusable pump sets are available with different size breast shields.

### 3.4 Disposable pump set – product parts (applied parts)

Symphony disposable pump sets are available with different size breast shields.

#### **i** Information

All pictures shown are for illustration purpose only. Actual product may vary due to product enhancement. Refer to the instructions for use provided with the Symphony pump sets for product description and instruction on handling.

See **chapter 11** for ordering information.

## 4. Cleaning

### CAUTION

- Unplug the breast pump before cleaning.
- These recommendations are not substitutes for official hospital procedures that may differ among institutions. Refer to hospital policy for cleaning of shared capital equipment.
- In the event of conflicting national cleaning and disinfection requirements, these shall prevail over the Medela recommendations.
- Surface disinfectant agents for medical devices will not cause mechanical damage to the Symphony breast pump.
- Follow the instructions provided by the supplier of the surface disinfectant agent.
- Provide personal protection for operators in accordance with the supplier's instructions and safety data sheets.
- Check the suitability of the surface disinfectant agents for medical devices by reference to the supplier's information and/or physical testing.

### 4.1 Breast pump – cleaning and disinfection

Clean and disinfect the breast pump as needed.

1. To clean the breast pump, wipe the pump housing with a clean, damp (not wet) cloth. Drinking-quality water and a basic washing-up liquid or non-abrasive detergent may be used.

2. To disinfect the pump, wipe the pump housing with a clean, damp (not wet) cloth that is impregnated with a surface disinfectant agent for medical devices. Medela recommends using an alcohol-based surface disinfectant agent.

### 4.2 Cleaning pump sets

Refer to the instructions for use provided with the Symphony pump sets for product description and instruction on handling.

## 5. Preparing for pumping

### WARNING

- Only use the power cord that comes with the Symphony breast pump.
- Make sure the voltage of the power cord is compatible with the power source.

### Information

#### Optional for pump with battery:

- The built-in battery must be charged for 12 hours on receipt of the device.
- The pump can be used during charging.
- If the pump is to remain unused for a long time, charge the battery every two months (for 12 hours each time).
- The optimum service life of the battery is achieved if it is not fully discharged on a regular basis or if the pump is charged while not in use


#### Pumping time (fully charged battery):

approx. 60 minutes

#### Charging time (for full charge):


12 hours

### 5.1 Installation of the program card

1. The program card must be inserted in the card compartment on the back of the breast pump.  
→ For more information on the use of the program card and the selection of the appropriate program, see **section 6.2**.


### 5.2 Installation of the card and cord protector (optional)


1. Insert the cord into the card and cord protector in a loop around the pin.  
→ Ensure that there is sufficient cord length to insert the power cord into the power socket in the back of the pump.

## Preparing for pumping (cont.)

16


2. Push the card and cord protector over the guide rail on the back until it stops and tighten the screw.


3. To remove, unscrew the card and cord protector and slide out of the guide rail.

### 5.3 Mains operation

1. Insert the power cord into the power socket on the back of the pump.

2. Plug the power cord into the mains socket.


### 5.4 Assembling the Symphony pump set

#### CAUTION

- Wash hands thoroughly (at least one minute) with soap and water before touching the breast, the clean pump and the clean pump set parts.
- Dry your hands with a fresh towel or a single-use paper towel.


## NOTICE

- Use Medela original accessories only.
- Refer to instructions for use provided with the Symphony pump set for description and instruction on handling the pump set.
- Check pump set parts for wear or damage before use. Replace damaged or worn parts.
- All parts must be completely dry before use.


1. Carefully assemble the Symphony pump set according to the instructions for use provided with the pump set.
2. Screw a bottle into the connector. Use a bottle according to the individual needs (various bottle sizes available).


## 5.5 Connecting the pump set to the pump housing

1. Open the lid of the breast pump: Press the oval button on the top of the Symphony pump and lift the handle.

2. Position the Symphony cap firmly on the breast pump membrane to form a seal.  
→ Make sure that the Symphony cap engages and lies flush to the surface.
3. Centre the tubing connector in the recess in the pump housing.

### **i** Information

Always close the lid when pumping. The lid keeps the membrane cap pushed down to prevent a loss of vacuum that could occur if the membrane cap came off.

4. Close the lid.

## 6. Expression

18

### 6.1 Selecting the correct Medela breast shield size

Make sure to use the correct breast shield size. It is an essential component of effective pumping, helping to optimise milk flow.

#### Understanding Medela's breast shield sizing

Pumping should not hurt. For maximum comfort and pumping efficiency, Medela offers five breast shield sizes.

This guide is a starting point to help determine the optimal size based on the nipple diameter.


##### • Step 1

Using a ruler or measuring tape, measure the diameter of the nipple at the base (across the middle) in millimeters (mm). Do not include the areola.

##### • Step 2

Based on the measurement, determine the Medela breast shield size.

Example: If the nipple size measures 16 mm in diameter, the recommended Medela breast shield size is 21 mm.


See **chapter 11** for ordering information.

## Optimising the breast shield size


- Centre the nipple and gently hold the breast shield against the breast.
- Adjust for Maximum Comfort Vacuum to achieve optimum suction level.
- Refer to the images below while pumping in the expression phase (after the stimulation phase).

- The nipple rubs along the side of the tunnel.
- Try a larger size.

- The nipple is centred and moves freely.

- The nipple and excessive areola are pulled into the tunnel.
- Try a smaller size.

## Reasons for trying a new size

- If the nipple rubs against the tunnel sides to the point of causing discomfort.
- If excessive areola is being pulled into the tunnel.
- If there is any redness.
- If the nipple or areola is turning white.
- If the mother feels unexpressed milk after pumping.
- As levels of tissue edema subside post birth.

If any of these conditions occur, consider trying a different size by following the measuring instructions above.

## 6.2 Selecting the right program

The Symphony breast pump is installed with the Symphony PLUS program card, which contains two pump programs.

The goal of the INITIATE program is to initiate milk production. The MAINTAIN program is designed to build and maintain milk production.

Until secretory activation (milk coming in) or a maximum of 5 days:

Use the INITIATE program directly after birth until the mother has reached secretory activation:

- **For fully pump-dependent mothers:**

This can be indicated by pumping at least 20 ml in total in each of her last three pumping sessions. If not all pumping volumes have been recorded, this can be indicated by other physiological markers of secretory activation, such as breast fullness.

- **When the mother is breastfeeding and pumping:**

This can be indicated by other physiological markers of secretory activation, such as breast fullness

OR

- Throughout the first five days. If milk production has not started after five days, change to the MAINTAIN program.

After secretory activation (milk coming in):

Use the MAINTAIN program after secretory activation:

- Once the mother has pumped at least 20 ml in total in each of her last three pumping sessions, OR
- Beginning on day six, whichever of the two options comes first.


For more detailed information on the program card refer to the card's instructions for use.


## 6.3 Operation: INITIATE program

### ⚠ CAUTION

- Wash hands thoroughly (at least one minute) with soap and water before touching the breast, the clean pump and the clean pump set parts. Dry your hands with a fresh towel or a single-use paper towel.

1. Press the on/off button  to switch on the breast pump.  
→ The display text shows "INITIATE press ".

2. Press the "Let-down" button  within ten seconds.  
→ The display text will change to "INITIATE running".

3. Set the maximum comfort vacuum (see **section 6.5**).  
→ The bullets on the display show the vacuum level.  
The more bullets, the stronger the vacuum.

4. The INITIATE program runs automatically for 15 minutes with stimulation, expression and pause phases.  
→ An acoustic signal indicates the end of the program. The display shows "Program complete". The pump switches off automatically.  
It is important to complete the entire INITIATE program.

### NOTICE

- Always disconnect the Symphony breast pump from the power source after expression.


## 6.4 Operation: MAINTAIN program

### ⚠ CAUTION


- Wash hands thoroughly (at least one minute) with soap and water before touching the breast, the clean pump and the clean pump set parts. Dry your hands with a fresh towel or a single-use paper towel.

### i Information


- If the pump has been running for 30 minutes without any manipulations (e.g. adjustment of the vacuum), the Symphony breast pump will automatically switch off.

1. Press the On/off button  to switch on the breast pump.  
→ The MAINTAIN program will automatically start with the stimulation phase after pressing the On/off button.


2. Set the maximum comfort vacuum (see **section 6.5**).  
→ The drops on the display show the vacuum level.  
The more drops, the stronger the vacuum.

3. The pump will automatically switch to the expression phase after two minutes.
4. If the milk starts to flow sooner, press the Let-down button  to change to the expression phase.

5. Set the maximum comfort vacuum (see **section 6.5**).  
→ The bars on the display show the vacuum level.  
The more bars, the stronger the vacuum.

6. Continue pumping until the breast feels well drained and soft all over (Medela recommends a pumping time of at least 15 minutes).
7. Switch off the pump by pressing the on/off button .

## **i** Information

- Press the “Let-down” button  to revert from the expression phase to the stimulation phase if necessary. If milk has not begun to flow after two successive stimulation phases, the mother should take a break and try again in 15 to 30 minutes. Also, she could try breast massage and then pump using the expression phase.

## NOTICE

- Always disconnect the Symphony breast pump from the power source after expression.

### 6.5 Set maximum comfort vacuum

Set the maximum comfort vacuum by rotating the vacuum adjustment knob.


1. First rotate to the right to increase the vacuum until there is a slight discomfort.
2. Then rotate to the left to reduce the vacuum until it is comfortable.

### 6.6 Expression by single pumping

#### i Information


Medela recommends double pumping.

- Double pumping saves time and increases the energy content of milk. Milk production may be increased and maintained long term.
- For single pumping, step **7** of **section 6.7** can be omitted.
- Do not kink the tubing during expression.


1. Prepare one bottle and pump set for expression.
2. Mount only one Symphony cap on the breast pump (see **section 5.5**).
3. Continue as described in **section 6.7 step 3** for double pumping.


## 6.7 Expression by double pumping


1. Prepare two bottles and pump sets for expression.
2. Mount the Symphony caps on the breast pump (see **section 5.5**).

3. Place the breast shield on the breast so that the nipple is properly centred in the tunnel.
4. Hold the breast shield onto the breast with the thumb and index finger. Support the breast with the palm of the hand.


5. Switch the breast pump on with .
6. Proceed with the program selection (see **sections 6.2, 6.3 and 6.4**).

7. Place the breast shield of the second bottle on the other breast so that the nipple is properly centred in the tunnel.
8. Operate the pump as appropriate for the selected program (see **sections 6.3 and 6.4**).

## Expression (cont.)

26

9. When expression is finished, place the bottles in the stand or bottle holder in order to prevent the bottles from tipping over.

10. Close the bottles with the lid.

## 7. Storing and thawing breast milk

### 7.1 Storage and thawing guidelines

Use of a medical grade refrigerator with temperature uniformity and constant temperature monitoring is recommended.

Refer to [www.medela.com](http://www.medela.com) for a recommendation on storing and thawing of breast milk.

You can contact your local customer service in order to get a printed version if requested.

## 8. Troubleshooting and error codes

### 8.1 Troubleshooting

Problem	Solution
<b>If the motor is not running</b>	<p>Check whether mains power/battery is available.</p> <p>Check that the device is switched on.</p> <p>Check whether program card is inserted correctly in the card compartment on the back of the breast pump.</p>
<b>If there is low or no suction</b>	<p>Assemble the pump set in accordance with the instructions in <b>section 5.4</b>.</p> <p>Check whether all connections to the pump set and breast pump are secure.</p> <p>Express in accordance with <b>chapter 6</b>. Make sure that the breast shield forms a complete seal around the breast.</p> <p>Check that the pump set components are clean and not damaged.</p> <p>Check that the pump set is correctly assembled according to <b>section 5.4</b> and <b>5.5</b>:</p> <ul style="list-style-type: none"><li>• Check that the ends of the tubing are firmly connected.</li><li>• Do not kink or trap the tubing during expression.</li><li>• All components must be completely dry prior to pumping.</li></ul>
<b>If there is an overflow</b>	<ul style="list-style-type: none"><li>• Place the pump set into the bottle holder.</li><li>• Turn off the pump and unplug the power cord from the socket.</li><li>• Use a damp (not wet) cloth to wipe the pump and membrane area.</li><li>• If there is milk, condensation or visible residues in the tubing:<ul style="list-style-type: none"><li>For disposable pump set: Dispose of the pump set and replace with a new one.</li><li>For reusable pump set: Dispose of the silicone tubing and replace with a new one. Reprocess the remaining parts according to <b>section 4.2</b>.</li></ul></li><li>• When components are completely dry, reassemble breast pump and accessories set.</li></ul>

# Troubleshooting and error codes (cont.)

## 8.2 Error codes

The following table describes error codes that may be displayed indicating an issue with the unit, potential causes, and steps to correct the issue.

Display	Potential issues	Actions
Invalid card	Program card is not a Symphony card Program was inserted incorrectly Program card is damaged and not working	<ul style="list-style-type: none"><li>• Check if the correct program card is being used</li><li>• Make sure that the card is positioned correctly</li><li>• Replace program card</li></ul>
No card	Program card is not inserted Program card is not completely inserted	<ul style="list-style-type: none"><li>• Insert program card</li><li>• Push card all the way in</li></ul>
Battery low	Battery needs charging soon, will beep every 20 seconds (15 minutes of pumping time remaining from when the indication started)	<ul style="list-style-type: none"><li>• Make sure that the power cord is plugged into the back of the pump</li></ul>
Motor overload	Motor is drawing too much current	<ul style="list-style-type: none"><li>• Contact your local Medela representative if this happens repeatedly</li></ul>
Error	Potential fault in electronics of gear control	<ul style="list-style-type: none"><li>• Contact your local Medela representative if this happens repeatedly</li></ul>
Power fault	Potential fault in electronics of power supply	<ul style="list-style-type: none"><li>• Contact your local Medela representative if this happens repeatedly</li></ul>

## 9. Warranty and maintenance

### Warranty

Three-year warranty for the unit from the date of purchase, excluding pump set and accessories. Six month warranty for the batteries. The manufacturer is not liable for any damage or consequential damage caused by incorrect operation, inappropriate usage or use by unauthorised persons. Routine checks, servicing and battery exchange may only be carried out by Medela-authorized personnel.

### Electrical safety

The engineering of the Medela Symphony breast pump has been verified by independent test institutes for compliance with the EN/IEC 60601-1 standard. Copies of such certificates may be obtained upon request. Because of the method of construction used, Medela does not expect electrical safety to be affected at any time during the life of the product – provided that it is used properly in accordance with the intended use as well as these instructions for use, and that the device is serviced and repaired at authorised service centres. There are no prescribed or recommended electrical safety checks or recurrent tests to be carried out.

Manufacturer's information provided in accordance with EN/IEC 62353:2014 "Medical electrical equipment – Recurrent test and test after repair of medical electrical equipment".

### Maintenance

Besides the cleaning procedure(s) described in **chapter 4** and battery charging for pumps with optional batteries as described in **chapter 5**, there is no further maintenance required. Besides the inspection procedures described in **chapter 1**, e.g. visual signs of damage on pump or power cord, pump dropped or dropped into water, there are no further inspections required. As pointed out in **chapter 1**, service and repair must be performed only by an authorised service agency (centre) in full accordance with the Symphony Service Manual using original Medela spare parts mentioned therein. The service manual, wiring diagrams, and descriptions are available upon request from Medela.

Medela deems that the final production testing replaces on-site testing of the device before putting into service.

The service life of the device is seven years, internal batteries life included.

The service life of the accessories is two years.

# 10. Disposal

30

## Disposal of electrical and electronic equipment

The Symphony breast pump is made of metals and plastics and should be disposed of in accordance with the European directives 2011/65/EU and 2012/19/EU. Additionally, local guidelines must also be observed. With the AC/DC versions, the electronic components and the rechargeable battery must be disposed of separately, in accordance with the local regulations. Please take care that you dispose of the Symphony and its accessories in accordance with your local disposal guidelines.

### User information for the disposal of electrical and electronic equipment

This symbol means that the electrical and electronic equipment must not be disposed of together with unsorted municipal waste. Correct disposal of this device protects and prevents possible damage to the environment or human health. For more information about the disposal contact the manufacturer, your local caregiver or healthcare provider. This symbol is only valid in the European Union. Please respect the relevant state laws and rules in your country for the disposal of electrical and electronic equipment.

# 11. Accessories list

## Symphony PLUS program card

Article number	Product
200.9186	Symphony PLUS program card English
200.9187	Symphony PLUS program card German
200.9188	Symphony PLUS program card French
200.9189	Symphony PLUS program card Dutch
200.9190	Symphony PLUS program card Spanish
200.9191	Symphony PLUS program card Italian
200.9192	Symphony PLUS program card Portuguese
200.9193	Symphony PLUS program card Swedish
200.9194	Symphony PLUS program card Norwegian
200.9195	Symphony PLUS program card Danish
200.9196	Symphony PLUS program card Finnish
200.9197	Symphony PLUS program card Polish
200.9198	Symphony PLUS program card Japanese

## Symphony Standard program card

Article number	Product
800.0540	Symphony Standard program card German
800.0541	Symphony Standard program card English
800.0542	Symphony Standard program card French
800.0543	Symphony Standard program card Italian
800.0544	Symphony Standard program card Dutch
800.0545	Symphony Standard program card Swedish
800.0546	Symphony Standard program card Danish
800.0547	Symphony Standard program card Finnish
800.0548	Symphony Standard program card Norwegian
800.0549	Symphony Standard program card Polish
800.0550	Symphony Standard program card Spanish
800.0551	Symphony Standard program card Portuguese
800.0552	Symphony Standard program card Greek
800.0553	Symphony Standard program card Hungarian
800.0554	Symphony Standard program card Czech

Articles may not be available in all countries. For further information about Medela products, visit [www.medela.com](http://www.medela.com) and look for the location finder.

## Accessories list (cont.)

32

### Disposable single-use pump sets

Article number	Product
008.0397	Sterile single-use Symphony pump set size M (24 mm)
008.0398	Sterile single-use Symphony pump set size L (27 mm)
008.0399	Sterile single-use Symphony pump set size XL (30 mm)
101033930	Symphony Pump Set with PersonalFit PLUS, sterile, S (21mm) (30 pcs)
101033931	Symphony Pump Set with PersonalFit PLUS, sterile, M (24mm) (30 pcs)
101033932	Symphony Pump Set with PersonalFit PLUS, sterile, L (27mm) (30 pcs)
101033933	Symphony Pump Set with PersonalFit PLUS, sterile, XL (30mm) (30 pcs)

### Disposable one-day pump sets

Article number	Product
008.0276	Ready-to-Use one-day Symphony pump set, size M (24 mm)
008.0277	Ready-to-Use one-day Symphony pump set, size L (27 mm)
008.0278	Ready-to-Use one-day Symphony pump set, size XL (30 mm)
101033934	Symphony Pump Set with PersonalFit PLUS, RtU, S (21mm) (30 pcs)
101033935	Symphony Pump Set with PersonalFit PLUS, RtU, M (24mm) (30 pcs)
101033936	Symphony Pump Set with PersonalFit PLUS, RtU, L (27mm) (30 pcs)
101033937	Symphony Pump Set with PersonalFit PLUS, RtU, XL (30mm) (30 pcs)

### Disposable single-use bottles

Article number	Product
200.2110	Colostrum container, set of 2 disposable bottles 35 ml with lids (50 pcs)
200.2957	Ready-to-Use single-use bottle, 80 ml, with lid (480 pcs)
200.2956	Ready-to-Use single-use bottle, 150 ml, with lid (320 pcs)
200.4802	Ready-to-Use single-use bottle, 250 ml, with lid (160 pcs)
008.0297	Sterile single-use bottle, 80 ml, with lid (120 pcs)
008.0304	Sterile single-use bottle, 150 ml, with lid (85 pcs)
008.0305	Sterile single-use bottle, 250 ml, with lid (60 pcs)

Articles may not be available in all countries. For further information about Medela products, visit [www.medela.com](http://www.medela.com) and look for the location finder.

## Reusable pump sets

Article number	Product
008.0486	Symphony pump set hospital multibox (S) (25 pcs)
008.0465	Symphony pump set hospital multibox (M) (25 pcs)

## Other accessories

Article number	Product
008.0467	PersonalFit PLUS breast shield S (21 mm) (25 pcs)
008.0468	PersonalFit PLUS breast shield M (24 mm) (25 pcs)
008.0469	PersonalFit PLUS breast shield L (27 mm) (25 pcs)
008.0470	PersonalFit PLUS breast shield XL (30 mm) (25 pcs)
008.0471	PersonalFit PLUS breast shield XXL (36 mm) (25 pcs)
008.0466	Connector kit (25 pcs)
008.0472	Membrane (25 pcs)
008.0473	Tubing (25 pcs)
008.0474	Symphony cap (25 pcs)
800.0625	Breast milk bottle 150 ml (25 pcs)
008.0475	Multi lid (25 pcs)
800.0627	Lid with opening (25 pcs)
800.0628	Disc (25 pcs)
800.0629	Cap (25 pcs)
600.0832	Mobile stand for Symphony complete
810.0552	Bottle stand for Symphony
600.0973	Hard case for Symphony with strap
600.7073	Hard case for Symphony without strap
600.0816	Symphony card and cord protector
928.1017	Car connection cable for Symphony
600.0387	Test vacuum meter

Articles may not be available in all countries. For further information about Medela products, visit [www.medela.com](http://www.medela.com) and look for the location finder.

# 12. Meaning of symbols

## Symbols used in these instructions

General safety alert symbol, points to information related to safety.

Identifies the manufacturer.

Defines a temperature range (e.g. for operation, transport or storage).

Defines a relative humidity range (e.g. for operation, transport or storage).

Defines an atmospheric pressure range (e.g. for operation, transport or storage).

Indicates compliance with the EU requirements concerning medical devices.

## Symbols on the display

Vacuum level indicator, stimulation phase – Standard program card and Symphony PLUS program card (MAINTAIN program).

Vacuum level indicator, expression phase – Standard program card and Symphony PLUS program card (MAINTAIN program).

Vacuum level indicator – Symphony PLUS program card (INITIATE program).

## Symbols on the device and packaging

On/off button

Let-down button

Read and follow the instructions for use.

Indicates compliance with the EU requirements concerning medical devices.

Identifies the manufacturer.

Do not dispose of electric/electronic devices together with unsorted municipal waste (dispose of the device in accordance with local regulations).



Indicates compliance with additional USA and Canada safety requirements for medical electrical equipment.

Authorized representative in the European Community / European Union.

Importer

	Indicates compliance with international requirements for protection from electric shock (Type BF applied parts).
	Indicates the serial number of the device.
	Indicates the part number of the device.
<b>IP21</b>	Indicates the degree of protection against ingress of foreign objects and moisture.
	Indicates the date of manufacturing.
	Indicates equipment meeting the safety requirements specified for Class II equipment according to IEC 61140.
	Indicates the fuse (to identify fuse boxes or their location).
	Indicates that the fuse box is safety tested.
	Indicates that the cover of the mains socket should be opened.
	Indicates a medical device.
	Indicates that the material is part of a recovery/recycling process.
	Indicates that the package is made of cardboard.
	Keep away from sunlight.
	Contains fragile goods. Handle with care.
	Defines a temperature range (e.g. for operation, transport or storage).
	Defines a relative humidity range (e.g. for operation, transport or storage).
	Defines an atmospheric pressure range (e.g. for operation, transport or storage).
	Keep away from rain. Keep in dry conditions.
	Indicates that the packaging contains products intended to come in contact with food according to regulation (EC) No 1935/2004.
	Indicates unique GSI Global Trade Item Numbers (GTIN).

## 13. Electromagnetic compatibility (EMC)

The Symphony breast pump is EMC-tested in conformity with the requirements of IEC 60601-1-2:2007 and 60601-1-2:2014 4th Edition according to clause 7 and 8.9.

### **WARNING**

- The Symphony breast pump should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the Symphony electric breast pump should be monitored to verify normal operation in the configuration in which it will be used.
- Wireless communications equipment such as wireless home network devices, mobile phones, cordless telephones and their base stations, and walkie-talkies can affect the Symphony electric breast pump and should be kept at a distance of at least 30 cm from the equipment.
- Do not use other accessories than those specified or sold by the manufacturer as replacement parts for internal components as it may result in increased emissions or decreased immunity of the Symphony device.

For technical details (English only) refer to section 15.1 in the appendix.

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

# 15. Appendix

## 15.1 EMC/technical description (English only)

The Symphony breast pump needs special precautions regarding EMC (Electromagnetic Compatibility) and needs to be installed and put into service according to the EMC information provided in these instructions for use.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

### NOTICE

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorientate or relocate the receiving antenna.
- Increase the distance between the equipment and receiver.
- Connect the equipment to a different outlet on the circuit than that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Electromagnetic Compatibility (EMC, IEC 60601-1-2:2007 3rd Edition and IEC 60601-1-2:2014 4th Edition, Table 1)

<b>Electromagnetic emissions</b>		
<p>The Symphony electric breast pump is intended for use in the electromagnetic environment specified below. The customer or the user of the Symphony electric breast pump should ensure that it is used in such an environment.</p>		
<b>Emission tests</b>	<b>Compliance</b>	<b>Electromagnetic environment – guidance</b>
<b>RF emissions CISPR 11</b>	Group 1	The Symphony electric breast pump uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
<b>RF emissions CISPR 11</b>	Class B	The Symphony electric breast pump is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
<b>Harmonic emissions IEC 61000-3-2</b>	Class A	
<b>Voltage fluctuations/flicker emissions IEC 61000-3-3</b>	Complies	

# Appendix (cont.)

## Electromagnetic Compatibility (EMC, IEC 60601-1-2:2007, Table 2)

<b>Electromagnetic immunity</b>			
The Symphony breast pump is intended for use in the electromagnetic environment specified below. The customer or the user of the Symphony electric breast pump should ensure that it is used in such an environment.			
<b>Immunity test</b>	<b>IEC 60601 test level</b>	<b>Compliance level</b>	<b>Electromagnetic environment – guidance</b>
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input / output lines	± 2 kV for power supply lines ± 1 kV for input / output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% $U_T$ (>95% dip in $U_T$ ) for 0.5 cycle 40% $U_T$ (60% dip in $U_T$ ) for 5 cycles 70% $U_T$ (30% dip in $U_T$ ) for 25 cycles <5% $U_T$ (>95% dip in $U_T$ ) for 5 s	<5% $U_T$ (>95% dip in $U_T$ ) for 0.5 cycle 40% $U_T$ (60% dip in $U_T$ ) for 5 cycles 70% $U_T$ (30% dip in $U_T$ ) for 25 cycles <5% $U_T$ (>95% dip in $U_T$ ) for 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Symphony electric breast pump requires continued operation during mains power interruptions, it is recommended that the Symphony electric breast pump be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
<b>Note:</b> $U_T$ is the a.c. mains voltage prior to application of the test level			

## Electromagnetic Compatibility (EMC, IEC 60601-1-2:2007, Table 4)

### Electromagnetic immunity

The Symphony breast pump is intended for use in the electromagnetic environment specified below. The customer or the user of the Symphony electric breast pump should ensure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	<p>Portable and mobile RF communications equipment should be used no closer to any part of the Symphony electric breast pump, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance  <math>d = 0.35\sqrt{P}</math>  <math>d = 1.2\sqrt{P}</math> 80 MHz to 800 MHz  <math>d = 2.3\sqrt{P}</math> 800 MHz to 2.5 GHz            where <math>P</math> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and <math>d</math> is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,<sup>a</sup> should be less than the compliance level in each frequency range.<sup>b</sup></p> <p>Interference may occur in the vicinity of equipment marked with the following symbol: </p>
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	

**Note 1:** At 80 MHz and 800 MHz, the higher frequency range applies.

**Note 2:** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

<sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Symphony electric breast pump is used exceeds the applicable RF compliance level above, the Symphony electric breast pump should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Symphony electric breast pump.

<sup>b</sup> Field strengths should be less than 10 V/m over the 150 kHz to 80 MHz frequency range.

Electromagnetic Compatibility (EMC, IEC 60601-1-2:2007, Table 6)

**Recommended separation distances between portable and mobile RF communications equipment and the Symphony electric breast pump**

The Symphony breast pump is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Symphony electric breast pump can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Symphony electric breast pump as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter <i>W</i>	Separation distance according to frequency of transmitter <i>M</i>		
	150 kHz to 80 MHz <i>d</i> = 0.35√ <i>P</i>	80 MHz to 800 MHz <i>d</i> = 1.2√ <i>P</i>	800 MHz to 2.5 GHz <i>d</i> = 2.3√ <i>P</i>
0.01	0.04	0.12	0.23
0.1	0.11	0.37	0.74
1	0.35	1.2	2.3
10	1.11	3.7	7.4
100	3.50	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance *d* in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where *P* is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

- Note 1:** At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.
- Note 2:** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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