Swing®

EN Instructions for use | ES Instrucciones de uso | FR Mode d'emploi

Important Safeguards

When using electrical products, especially when children are present, basic safety precautions should always be followed.

READ ALL INSTRUCTIONS BEFORE USING THIS PRODUCT.

DANGER: To reduce the risk of electrocution:

- Always unplug electrical product immediately after use.
- Do not use while bathing.
- Do not place or store product where it can fall or be pulled into a tub or sink.
- Do not place or drop into water or other liquid.
- Do not reach for a product if it has fallen into water. Immediately unplug from electrical outlet.

WARNING: To avoid fire, electrocution, or serious burns:

- Do not leave product unattended when plugged into an electrical outlet.
- Close supervision is necessary when this product is used near children or persons with disabilities.
- Use the product only for its intended use as described in this manual. Do not use attachments not recommended by the manufacturer.
- Never operate this product if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into water.
- Keep cord away from heated surface.
- Never use while sleeping or when drowsy.
- Never drop or insert any object into any opening or tubing.
- Do not use outdoors, or operate where aerosol (spray) products are being used where oxygen is being administered.
- Always inspect power adaptor and battery pack wires prior to use for damage or exposed wire. If damage is found, immediately discontinue use of power adaptor or battery pack and call Medela Customer Service at 1-800-435-8316.

The warning symbol identifies all instructions that are important to safety. Failure to observe these instructions can lead to injury or damage to the breast pump. When used in conjunction with the following words, the warning symbols stand for:

- WARNING** Can lead to serious injury or death.
- CAUTION** Can lead to minor injury.
- NOTE** Can lead to material damage.
- INFO** Useful or important information that is not related to safety.

WARNING: To avoid health risk and reduce the risk of injury:

- This product is intended for use by a single user only and should not be shared between users.
- Do not drive while pumping.
- Never use while pregnant, as pumping can induce labor.
- Clean and sanitize all parts that come in contact with your breast and breast milk prior to first use.
- Wash all parts that come in contact with your breast and breast milk after every use.
- Inspect all appropriate pump components before each use.
- Do NOT continue pumping for more than 2 consecutive pumping sessions if no results are achieved.
- Do not thaw frozen breast milk in a microwave or in a pan of boiling water.
- If tube becomes moldy, discontinue use and replace tube.
- If breast milk backup occurs you must clean tube before your next pumping session.

CAUTION: Can lead to minor injury:

- Do NOT wrap cord around adaptor body.
- Use only the power adaptor that comes with the Swing Breast Pump.
- Pump only with the breast pump in an upright position.
- Do NOT shorten tubing.
- Make sure the voltage of the power adaptor is compatible with the power source.
- Plug the power adaptor into the breast pump first and then into the wall socket.
- Do NOT use antibacterial or abrasive cleaners/detergents when cleaning breast pump or breast pump parts.
- Never put breast pump motor unit in water or a sterilizer, as you can cause permanent damage to the breast pump.
- Do not attempt to remove the breast shield from your breast while pumping. Turn the breast pump off and break the seal between your breast and breast shield with your finger, then remove breast shield from your breast.
- If pumping is uncomfortable or causing pain, turn the unit off, break the seal between the breast and the breast shield with your finger and remove the breast shield from your breast.
- Plastic bottles and component parts become brittle when frozen and may break when dropped.
- Bottles and component parts may become damaged if mishandled, e.g. dropped, over-tightened or knocked over.
- Take appropriate care in handling bottles and components.
- Do not use the breast milk if bottles or components become damaged.

SAVE THESE INSTRUCTIONS

Questions? Visit www.medela.com or call us at 1-800-435-8316.

1. Cleaning	8-13
1.1 Before using for the first time	8
1.2 After each use	8-9
1.3 Sanitize daily.....	10-11
1.4 Tubing and motor unit cleaning	12-13
2. Assembly	14-21
2.1 Assembly of pump kit	14-17
2.2 PersonalFit™ Sizing	18
2.3 Powering your pump	19-21
2.3.1 Power adaptor care	19
2.3.2 Battery compartment.....	20-21
3. Pumping	22-26
3.1 Device operation.....	22-25
3.2 Mobility.....	26
4. Storing, Preparing & Feeding Breast Milk	27
4.1 Storing Breast Milk	27
4.2 Preparing Breast Milk.....	27
4.3 Feeding Breast Milk	27
5. Troubleshooting	28-29
6. Supplemental Information	30
7. EMC / Technical Description	31-37
8. Technical Specifications / Disposal Guideline	38
9. Warranty	39

Take a moment to read through this entire instruction manual before using this product for the first time.

PLEASE SAVE THESE INSTRUCTIONS.

It is best to wait until your breastfeeding routine is established (approximately 4 weeks) before expressing breast milk, unless otherwise advised by your healthcare professional.

Glossary of Terms

2-Phase Expression® technology – research based technology that mimics a baby’s natural nursing rhythm.

Expression Phase – slower sucking/pumping rhythm for gentle and efficient milk removal as quickly as possible.

Maximum Comfort Vacuum™ – highest vacuum setting where a mother still feels comfortable during pumping. It is different for every mother.

Stimulation Phase – fast sucking/pumping rhythm to stimulate the milk ejection reflex and to start the milk flowing.

Intended Use

For lactating women to express and collect milk from their breasts to complement breastfeeding.

If you have medical reasons or other needs for exclusively pumping, it is recommended that you use a hospital grade breast pump such as our Symphony Breast Pump. To learn more, visit us at www.medelabreastfeedingus.com.

Product Description

Swing® is a personal use electric breast pump that includes 2-Phase Expression technology and is capable of single pumping.

Product Description

1x PersonalFit™ Breast Shield (24mm)
(for other breast shield sizes see **Section 2.2**, page 18)
87073 (24 mm)

1x Connector
8007196

1x Valve
87089

3x Membranes
87088

2x Bottles
6007135

1x Bottle Stand
8100462

1x Neck Strap
8000734


1x Swing Storage Bag
8000838

1x Swing Tubing
8007215


Overview of LEDs

lights up	On
does not light up	Off
even, regular flashing	Stimulation phase
light constantly on	Expression phase
rapid flashing	Min./Max. vacuum

1x Motor Unit
0300011

1x Backplate

2x Lids
8107183


1x Power Adaptor
2011059

1. Cleaning

1.1 Before using for the first time

Clean and sanitize the parts below before using your pump for the first time. Follow the cleaning and sanitizing instructions in **Sections 1.2 (start at step 4)** and **1.3** that follow.

1.2 After each use


Clean all parts that come in contact with your breast and breast milk. See image in **Section 1.1**.

Separate breast shield from connector.

Twist and pull the yellow valve off the connector.

Remove the white membrane from the yellow valve.

Supplies Needed:

- Mild Dish Soap
- Clean Dish Towel or Soft Brush
- Clean Sink or Bowl
- Quality Drinking Water

Washing Options:

Washing in sink:

Rinse in cool water all separated parts that came in contact with breast and breast milk in order to remove breast milk residue.

- Soak all separated parts in warm soapy water for 5 minutes.
- Clean with a clean dish-cloth or soft brush.
- Rinse all separated parts with clear water.
- Allow all parts to air dry in a clean area.
- Store dry parts when not in use.

or

- You may wash your kit parts by following the instructions on the Quick Clean™ Breast Milk Removal Soap bottle.

Washing in dishwasher:

- Wash all separated parts on top rack of dishwasher.
- Allow all breast pump parts to air dry in a clean area.
- Store dry parts when not in use.

Caution

- Only use drinking-quality tap or bottled water for cleaning.
- Take apart and wash all parts that come in contact with the breast and breast milk immediately after use to avoid drying of milk residue and to prevent growth of bacteria.

1. Cleaning (cont.)

Note

- **Wash** breast pump kit parts after each use.
- **Sanitize** breast pump kit parts once daily.
- The pump kit requires no maintenance.
- Take care not to damage parts of the pump kit during cleaning.
- If the individual components of the pump kit are cleaned in a dishwasher parts may be discolored by food pigments. This has no impact on the function.
- Store the pump kit in a clean bag/container until next use.

1.3 Sanitize daily


Wash hands thoroughly.

Separate all parts that come in contact with your breast and breast milk.

- Fill a pot with enough water to cover all parts.
- Bring water to a boil.
- Place parts in boiling water for 10 minutes.

Allow water to cool and gently remove parts from water with tongs.

- Place parts on a clean surface and/or towel.
- Allow all parts to air dry.
- Store dry parts when not in use. Do NOT store wet or damp parts.

or

You may sanitize your kit parts by following the instructions on Medela's Quick Clean™ Micro-Steam™ bag. (Do not sanitize tubing in Micro-Steam bag.)

1. Cleaning (cont.)

Caution

- If tubing becomes moldy, discontinue use and replace tubing. To find replacement parts visit www.ShopMedela.com.
- Do not wash power adaptor.

Note

- Tubing should only be washed if dirty or milk is present. Follow instructions in **Section 1.4**.
- Tubing should be washed if condensation due to breast milk is present.
- Cleaning tubing is not necessary if condensation present is from previous washings or atmospheric conditions.

1.4 Tubing and Motor Unit Cleaning


Turn off breast pump.

Unplug breast pump from power source.

Remove tube from motor unit.

Remove tube from connector.

Rinse tube in cool water to remove breast milk.

- Wash tube in warm soapy water.
- Rinse tubing with clear water.
- Hang to air dry.

Wipe off with a clean, damp cloth.

2. Assembly


Parts you will need for this step:

2.1 Assembly of pump kit

Wash hands before touching breast pump, kit and breasts and avoid touching the inside of containers or lids.

Snap membrane onto valve until membranes lies completely flat on valve.

Push yellow valve securely onto inner ring of connector.

Screw bottle onto connector.

Insert PersonalFit™ breast shield into top opening on connector.

Insert hard yellow adaptor end of tube into connector.

Continued next page.

! Note

- Inspect valves for cracks or chips and inspect membranes for holes, tears or warping. If you notice any damage to these parts, replace them immediately.

i Info

- Carry out all the steps with care and assemble the pump kit correctly for optimum performance.

⚠ Caution

- Wash hands thoroughly with soap and water before touching the breast and pump parts.

2. Assembly (cont.)

Warning

- Use only the power adaptor supplied with your Swing Breast Pump.
- Always inspect the power adaptor prior to use for damaged, kinked, twisted, frayed or exposed wires. If damaged, immediately discontinue use and call Medela Customer Service at 1-800-435-8316.

Twists, kinks or bends

2.1 Assembly of pump kit (cont.)

Insert the triangle end of tube into triangle port on motor unit.

Plug the metal end of the power cord into the power adaptor port.

Plug power adaptor into a standard electrical outlet.

Caution

- While some discomfort may be felt when first using a breast pump, using a breast pump should not cause pain.
- If you are experiencing discomfort at the base of the nipple due to rubbing of your breast tissue with the breast shield tunnel, use of a lubricant such as Tender Care™ Lanolin may be beneficial.
- For assistance with correct breast shield sizing and comfort please visit www.medelabreastshields.com or see a lactation consultant/breastfeeding specialist.

Note

- See **Section 2.3.2** Battery Compartment for instructions on the use of the Battery Compartment.
- See **Section 2.3.3** Traveling Outside the U.S. for information on using your device internationally.

2. Assembly (cont.)

2.2 PersonalFit™ Sizing

1. Medela's PersonalFit sizing

Pumping should not hurt. For maximum comfort and pumping efficiency, Medela offers five breast shield sizes.

21 mm

24 mm

27 mm


30 mm

36 mm


Visit MedelaBreastShields.com to view the full guide to help you determine your optimal size based on your nipple diameter.

2. Test your breast shield size


1. Start with the 24 mm that came with your pump, or the size determined by measuring.
2. Center nipple and gently hold breast shield against your breast.
3. Adjust for Maximum Comfort Vacuum™, to achieve optimum suction level.
4. Refer to images while pumping in expression.

- Nipple rubs along side of tunnel
- Try a larger size

- Nipple is centered and moves freely

- Nipple and excessive areola are pulled into the tunnel
- Try a smaller size

3. Should you try a new size?

- Does your nipple rub sides of tunnel, to the point of causing discomfort?
- Do you see excessive areola being pulled into tunnel?
- Do you see any redness?
- Is your nipple or areola turning white?
- Do you feel unexpressed milk after pumping?

If you answered **“YES”** to any of these questions, consider trying a new size by following the measuring instructions above.

If you are still unsure if you selected the correct size, see a lactation consultant, breastfeeding specialist, or visit MedelaBreastShields.com for assistance in choosing the right size breast shield.


Caution

- While some discomfort may be felt when first using a breast pump, using a breast pump should not cause pain. If you are unsure about breast shield sizing, please contact a healthcare professional or breastfeeding specialist who can help you get a proper fit.

2.3 Powering your pump

2.3.1 Power Adaptor Care

To care for your adaptor be sure to follow these storage directions.

Do not wrap cord of power adaptor around the plug body.

Correct!

Do not unplug power adaptor by pulling on the cord.

Correct!

2. Assembly (cont.)

Note

- Remove the batteries if they are not going to be used for an extended period of time.
- If breast pump does not turn on, make sure the batteries have been inserted properly. If this does not solve the problem refer to the troubleshooting section of this instruction booklet or contact Medela Customer Service at 1-800-435-8316.
- Each set of batteries will provide approximately 2 hours of pumping time.
- Do not mix old and new batteries.
- Do not mix alkaline, standard (carbon-zinc) or rechargeable batteries.

2.3 Powering Options (cont.)

2.3.2 Battery Compartment

Use 4 AA batteries.

Replace battery cover.

Replace all batteries at the same time.

Remove battery cover by pushing the button toward the clip.

Place the NEGATIVE (-) end of the battery against the front of the spring by pushing the battery in and down in one motion until the battery snaps into place. The POSITIVE (+) end of the battery should be near the positive sign on the battery pack.

Note

- While travelling internationally, we recommend purchasing a Universal Power Plug Adapter (not included). Please consult with the country you are visiting to find out what type of adaptor will work best.

For information about traveling with your breast pump and expressed breast milk, visit www.TSA.gov.

3. Pumping

Caution

- Always inspect breast shield, connector, white membrane, yellow valve, and tube prior to use for cleanliness. If dirty, see **Section 1.2**.
- If you experience pain while pumping seek advice from your healthcare professional or breastfeeding specialist. See **Section 6** for more information.
- Do not try and express with vacuum that is too high and uncomfortable (painful). The pain, along with potential breast and nipple trauma, may decrease milk output.
- Make sure tube is not kinked or pinched while pumping.
- Contact your healthcare professional or breastfeeding specialist if you can express only minimal or no milk or if expression is painful.

Info

- Research has shown that you will pump more efficiently – get more milk in less time – when pumping at Maximum Comfort Vacuum during the Expression Phase.
- You should reassess your Maximum Comfort Vacuum throughout your pumping experience because it can change during the different stages of lactation.
- Do not tilt or overfill containers when pumping.
- Do not kink tubing.

3.1 Device Operation

1 Wash hands before touching breast pump, kit and breasts and avoid touching the inside of containers or lids.

2 Center the assembled breast pump and breast shield over your nipple.

3 Hold the breast shield onto your breast with your thumb and index finger. Support your breast with the palm of your hand.

4 Switch on the breast pump with . Your pumping session will start in the Stimulation Phase.

5 After two minutes, the pump will automatically change to the Expression Phase.

6 The vacuum can be adjusted in every phase. Comfort vacuum: Use to increase the vacuum until it is slightly uncomfortable. Then press to dial it back minimally. **Note:** The LED light will flash rapidly when you have reached minimum and maximum pumping settings.

Overview of LEDs

lights up	On
does not light up	Off
even, regular flashing	Stimulation phase
light constantly on	Expression phase
rapid flashing	Min./Max. vacuum

Continued next page.

3. Pumping (cont.)

Note

- Always disconnect Swing® from the power source after expression.
- Only fill the bottle to the 150 mL mark.
- Swing switches off automatically 30 minutes after last setting adjustment.

Info

- **How often should you pump?**
A breast pump is a replacement for when you are separated from your baby. It is important to pump when the baby would be breastfeeding, which usually means a working mother pumps 2-3 times during an 8 hour working day.
- **How long should your pumping session last?**
Pumping times can vary from mother to mother sometimes 15 minutes, sometimes up to 30 minutes.
- **How should your breasts feel after pumping?**
After pumping, your breasts should feel soft and there should be no firm areas. Before pumping, your breasts will have a firm, heavy feeling. After pumping make sure your breasts are not firm or have lumpy areas. This could indicate that the breast is not draining all over. See **Section 6** for more information.

3.1 Device Operation (cont.)

7 Switch off the breast pump with .

10 Disassemble and clean per Cleaning Instructions in **Section 1**.

8 Close the bottle with a lid. Follow the instructions in **Section 4** “Storing and thawing breast milk”. Use bottle stand to prevent the bottle from tipping over.


9 Unplug breast pump from power source.


3. Pumping (cont.)

Swing® is small and lightweight, and comes with everything you need.

3.2 Mobility

Use on any tabletop or flat surface.

Clip onto your belt or pocket.

Hang around your neck, over your shoulder, or from a hook.

4. Storing, Preparing and Feeding Breast Milk

4.1 Storing Breast Milk

Warning

- Do not thaw frozen breast milk in a microwave or in a pan of boiling water.
- Do not microwave breast milk. Microwaving can cause severe burns to baby's mouth from hot spots that develop in the milk during microwaving. (Microwaving can also change the composition of breast milk.)
- When freezing, do not fill container more than 75% full to allow space for expansion. Label container with date of pumping.

Freshly Expressed Breast Milk Storage Guidelines (For Healthy Term Babies)

Room Temperature	Cooler with Ice Pack	Refrigerator	Freezer
4–6 hours at 66–78 °F (19–26 °C)	24 hours at 59 °F (15 °C)	3–8 days at 39 °F or lower (4 °C)	6–12 months 0–4 °F (-18–20 °C)

References: www.BreastmilkGuidelines.com

4.2 Preparing Breast Milk

- Thaw breast milk overnight in the refrigerator. Thawed breast milk is safe in the refrigerator for 24 hours. Do not refreeze thawed breast milk.
- Quickly thaw breast milk by holding the bottle under warm running water.
- Place the sealed bottle in a bowl of warm water for 20 minutes to bring it to body temperature.
- If you are adding expressed breast milk to a container of already frozen breast milk, make sure to add a lesser amount than the already frozen amount.

4.3 Feeding Breast Milk

It is recommended that breastfeeding is well established prior to bottle feeding your baby.

- Always inspect the bottle, nipple and other components immediately before and after each use. If nipple appears cracked or torn, discontinue use immediately.
- To prevent possible choking hazard, test strength of nipple by pulling on bulb portion of the nipple.
- Do not attempt to enlarge the nipple hole.
- Infants must not be bottle fed without adult supervision.
- Nipple should not be used as a pacifier.


For additional Breast milk collection and storage information, please visit www.BreastmilkGuidelines.com

Questions? Visit www.medela.com or call us at 1-800-435-8316.


5. Troubleshooting

<p>Low or No Suction</p>	<p>Valve and Membrane Check</p> <ul style="list-style-type: none"> Remove and separate valve from membrane. Inspect the valve for chips or cracks. Inspect the membrane for holes or tears. If flawed or damaged replace valve and/or membrane prior to pumping. Wash valve and membrane per Cleaning section. Reassemble. Make sure the membrane lays flat on the valve, and if it does not, replace the valve and/or membrane prior to pumping. <p>Breast Shield and Connector Check</p> <ul style="list-style-type: none"> Inspect breast shield and connector to make sure there are no cracks and they securely attach to each other. Make sure the breast shield forms a complete seal around the breast. <p>Connection Point Check</p> <p>Inspect following connection points to ensure attachment is secure</p> <ul style="list-style-type: none"> Tubing adapter to connector. Tubing end to motor unit. Check that the tubing is not kinked or pinched. Other: Perform steps of troubleshooting Section "milk in motor". Test vacuum. If suction is not improved, contact Medela Customer Service at 1-800-435-8316.
<p>Power source</p>	<ul style="list-style-type: none"> Make sure pump is plugged in correctly. Pump will NOT turn on if tubing is plugged into the AC adapter port. Make sure the "on" button is being held between 1 and 3 seconds. Do not repeatedly hit the on/off button, or hold the on/off button down longer than 3 seconds. If using batteries, make sure transformer is not plugged into pump. If using batteries, make sure batteries are installed correctly inside pump. If pump is still not working with the power adaptor or the battery compartment, contact Medela Customer Service at 1-800-435-8316.
<p>Condensation or milk in tubing</p>	<ul style="list-style-type: none"> Refer to "Tubing and motor unit cleaning" section 1.4.
<p>No "let-down" or breast milk expression</p>	<ul style="list-style-type: none"> Ensure your breast pump is assembled correctly and that there is suction. Relax and take a 10-15 minute break if let down is not achieved after 2 consecutive pumping sessions. Consult with your healthcare provider or breastfeeding specialist if expression does not occur.

Milk in Motor

Turn off the pump and unplug from power source.

- Open the battery cover, remove the foam piece. Wash the foam piece with clear, lukewarm water and squeeze out well. Allow to dry thoroughly.
- Pour about 2 oz of warm soapy water into a container.
- Connect the tubing and the power cable onto the pump.
- Hold the device with the control button upwards over a basin and switch on. Press  and ensure that the lowest vacuum level is set by pressing  twice.
- Hold the end of the tubing that is not connected to the pump in the container with the soapy water.
- The water is pumped through the tubing into the device. Allow it to pump until all the water has been flushed through.
- Repeat procedure with 1 oz of clear, warm water.
- Then leave the device switched on for 3 minutes, until the pump is completely dry.
- As soon as the expanded-plastic liner is dry, put it in place again and close the battery cover.
- If necessary, clean the pump set according to cleaning instructions "Prior to first use and once daily".

Note

- This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:
 - Reorient or relocate the receiving antenna.
 - Increase the separation between the equipment and receiver.
 - Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
 - Consult the dealer or an experienced radio/TV technician for help.

6. Supplemental Information

Supplemental Information

Following are some common breastfeeding related conditions. If you experience any of these symptoms, contact a healthcare professional or breastfeeding specialist.

	Symptom	Potential Causes
Engorgement	Breasts are hard and uncomfortable, possibly with reddened areas.	Milk not draining properly from breast. Milk can drain back into the tissue causing areas to swell and become tender.
Blocked/Plugged Ducts	Area of breast will look reddened and may be tender to touch.	Milk not being drained from a specific duct. The area becomes 'clogged' up and milk is then prevented from flowing.
Mastitis	It often follows on from engorgement but can just occur suddenly. The initial symptoms are similar to flu – fatigue, headache and muscle aches, fever and localized breast tenderness. It can be in one or even both breasts and needs immediate treatment.	A bacterial infection in the breast tissue which is often accompanied by cracked nipples.

A Healthcare Professional or Breastfeeding Specialist can give you guidance in breastfeeding your baby and using a breast pump. Their guidance is valuable in terms of the lifetime of health benefits for both you and your baby. A listing of Lactation Professionals are available at www.medelalocator.com or call 1-800-TELL-YOU (1-800-835-5968).

If you have not resolved the problem with your breast pump or you have further questions, please contact Medela Customer Service at 1-800-435-8316 or email customer.service@medela.com.

7. EMC / Technical Description

The Swing breast pump needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in the instructions for use. Portable and mobile RF communications can affect the breast pump.

Note

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

7. EMC/Technical Description (cont.)

Guidance and manufacturer's declaration – electromagnetic emissions

This breast pump is intended for use in the electromagnetic environment specified below. The customer or the user of the breast pump should assure that it is used in such an environment.

Electromagnetic Compatibility (EMC, IEC 60601-1-2:2014 4th Ed.)

Emission tests	Compliance	Electromagnetic environment – guidance
RF Emissions CISPR 11	Group 1	The breast pump uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The breast pump is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not Applicable	

Warning

Warning – This Swing breast pump should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, this breast pump should be observed to verify normal operation in the configuration in which it will be used.

Guidance and manufacturer's declaration – electromagnetic immunity

This breast pump is intended for use in the electromagnetic environment specified below. The customer or the user of the breast pump should assure that it is used in such an environment. The Sonata Breast Pump has no essential performance but was tested for immunity to electromagnetic disturbances and passed using the following criteria:

1. No visible change in the operation of the breast pump.
2. The breast pump changes settings but returns automatically to previous settings.
3. The breast pump changes settings but can return to previous settings by intervention of the user.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8.kV, ± 15 kV air	± 8 kV contact ± 2 kV, ± 4 kV, ± 8.kV, ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV 100 kHz repetition frequency	± 2 kV 100 kHz repetition frequency	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 0,5 kV , ± 1 kV Line-to-line	± 0,5 kV , ± 1 kV Line-to-line	Mains power quality should be that of a typical commercial or hospital environment.

7. EMC/Technical Description (cont.)

Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0 % UT ; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315°	0 % UT ; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315°	Mains power quality should be that of a typical commercial or hospital environment. If the user of the breast pump requires continued operation during power mains interruptions, it is recommended that the breast pump be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m, 50 or 60 Hz	30 A/m, 50/60 Hz	It may be necessary to position the breast pump further from sources of power frequency magnetic fields or to install magnetic shielding. The power frequency magnetic field should be measured in the intended installation location to assure that it is sufficiently low.

NOTE U_T is the a.c. mains voltage prior to application of the test level.

Guidance and manufacturer's declaration – electromagnetic immunity

This breast pump is intended for use in the electromagnetic environment specified below. The customer or the user of the breast pump should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 V 0,15 MHz – 80 MHz 6 V in ISM and amateur radio bands between 0,15 MHz and 80 MHz 80% AM at 1 kHz	6V	Portable and mobile RF communications equipment should be used no closer to any part of the breast pump, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommend separation distance $d = 1,2\sqrt{P}$ $d = 1,2\sqrt{P}$ 80 MHz - 800 MHz $d = 2,3\sqrt{P}$ 800 MHz - 2.7 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).
Radiated RF IEC 61000-4-3	10 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz	10 V/m (minimum)	Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol. 

7. EMC/Technical Description (cont.)

Note 1 At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the breast pump is used exceeds the applicable RF compliance level above, the breast pump should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the breast pump.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the breast pump

This Swing electric breast pump is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the electric breast pump can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the electric breast pump as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.2\sqrt{P}$	80 MHz to 800 MHz $d = 1.2\sqrt{P}$	800 MHz to 2,7 GHz $d = 2.3\sqrt{P}$
0,01	0,06	0,04	0,07
0,1	0,18	0,11	0,22
1	0,58	0,35	0,22
10	1,84	1,11	2,21
100	5,83	3,50	7,00

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the breast pump including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

8. Technical Specifications / Disposal Guideline

Technical specifications

vacuum (approx.)
-15....-295 mmHg
-2....-39 kPa
45 ... 120 cpm

In	Out
100-240V - 50/60Hz 230-240V - 50Hz AUS/NZ only	4.8V --- 0.6/1.2A

4 x 1.5V, Alkaline
AA, Mignon, LR6

121 x 121 x 57 mm

300 g 
200 g 

Operation
+5 °C to +40 °C

Transport / Storage
-25 °C to +70 °C

Operation
15% to 93% RH

Transport / Storage
15% to 93% RH

Operation / Transport / Storage
70 to 106 kPa

9. Warranty

Caution

This is a single user product. Use by more than one person may present a health risk and performance questions and voids the warranty.

This product is warranted by Medela LLC. to the original retail purchaser to be free from defects in material and workmanship for the period of one year for pump mechanism (90 days for parts and accessories) from the date of purchase. Warranty can only be claimed in the country of purchase. In the event of a defect, Medela will repair or, at Medela's option, replace this product, without charge for such replacement, parts or labor. Purchaser shall bear all expense for returning this product to Medela. This warranty does not apply to any product used commercially or which has been subjected to misuse, abuse or alteration.

ANY AND ALL IMPLIED WARRANTIES, INCLUDING THE WARRANTY OF MERCHANTABILITY, ARE LIMITED TO A DURATION OF 90 DAYS FROM DATE OF PURCHASE. SOME STATES DO NOT ALLOW LIMITATIONS ON HOW LONG AN IMPLIED WARRANTY LASTS, SO THE ABOVE LIMITATIONS MAY NOT APPLY TO YOU. THIS WARRANTY GIVES YOU SPECIFIC LEGAL RIGHTS AND YOU MAY ALSO HAVE OTHER RIGHTS WHICH VARY FROM STATE TO STATE.

Before you make a claim under this warranty, it may save time and expense to call Medela Customer Service (toll free) at 1-800-435-8316. You may also call this number for additional information concerning this warranty. If you wish to make a claim under this warranty, you must return this product to Medela with a return authorization number received from Medela Customer Service, prepaid, together with your dated bill of sale or other proof of purchase and a brief statement of the problem to the following address:

Medela LLC – Returns, Door 4501
1101 Corporate Dr.
McHenry, IL 60050
ATTENTION: RETURNS

Medela Canada, Inc.,
4160 Sladeview Crescent, Unit #8,
Mississauga, Ontario,
Canada L5L 0A1
ATTENTION: RETURNS

Call first for authorization number. Returns not accepted without an authorization number.

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Printed in the USA. Medela wordmark and logo, 2-Phase Expression, Invented by Medela 2-Phase Expression Technology and Swing are registered in the U.S. Patent and Trademark Office. Quick Clean, Micro-Steam, Maximum Comfort Vacuum and PersonalFit are trademarks of Medela.
Medela AG
Distributed and Assembled by Medela LLC
Distribuido y armado por Medela LLC
Distribu  et assembl  par Medela LLC
Distributed in Canada by Medela Canada, Inc.