Vario

SUCTION PUMP

Vario 18 AC/DC

- EN Instructions for use
- DE Gebrauchsanweisung
- FR Mode d'emploi
- IT Istruzioni per l'uso
- NL Gebruiksaanwijzing
- SV Bruksanvisning
- FI Käyttöohjeet
- NO Bruksanvisning
- ES Instrucciones de uso
- PT Instruções de utilização

Table of contents

- Congratulations 4**
- 1 Warnings and safety instructions 5**
 - Warnings and safety instructions for healthcare professionals 6
 - Safety instructions 7
 - Warnings and safety instructions for patients and caregivers 8
- 2 Power supply and battery operation 10**
 - Battery low (only AC/DC versions) 10
 - Battery empty (only AC/DC versions) 10
- 3 Description 10**
 - Introduction 10
 - Intended use 10
 - Contraindications 11
 - Intended user 11
 - Intended Patient Population 11
 - Important note 11
 - Expected service life 11
- 4 Overview 12**
 - Definition of vacuum 12
 - Versions of the pump 12
 - Front of the pump 12
 - Bottom of device 13
 - Operating elements and battery status 13
 - Acoustic indicator 13
 - On/off switch AC versions 13
 - On/off switch AC/DC versions 13
 - On/off switch c/i versions 13
- 5 Installation 14**
 - Check initial delivery 14
- 6 Preparation for use 15**
 - Checks before use 15
 - Assembly of the basic configuration 15
 - Assembly of the REUSABLE collection system 16
 - Assembly of the DISPOSABLE collection system 17
- 7 Operating instructions 18**
 - Connect Vario to mains power 18
 - Check maximal vacuum for functional check 18
 - Changing vacuum level 19
 - Placing out of operation after use 19
- 8 Changing reusable jar 20**

| | | |
|-----------|--|------------|
| 9 | Changing disposable liner | 21 |
| 10 | Troubleshooting | 22 |
| | Motor not running..... | 22 |
| | Insufficient vacuum..... | 22 |
| | No LED lit (AC/DC version only)..... | 22 |
| | LED lit up (AC/DC version only)..... | 22 |
| | LED blinks slowly (AC/DC version only)..... | 22 |
| | LED blinks rapidly (AC/DC version only)..... | 22 |
| 11 | Replacing defective fuse | 23 |
| 12 | Setting interval times | 24 |
| 13 | Replacing the QuatroFlex | 25 |
| 14 | Replacing the rechargeable battery | 26 |
| 15 | Battery data and test | 27 |
| | Battery low..... | 27 |
| | Battery empty..... | 27 |
| | Duration of pumping..... | 27 |
| | Charging time (charged up to 100%)..... | 27 |
| | Testing the rechargeable battery..... | 27 |
| 16 | Battery calibration | 28 |
| | Calibrating NiMH rechargeable battery..... | 28 |
| | Storage of the rechargeable battery..... | 28 |
| 17 | Filter test | 29 |
| 18 | Cleaning guidelines | 30 |
| | General notes..... | 30 |
| | Medela recommended surface cleaning agents for pump housing..... | 30 |
| | Water..... | 30 |
| | Cleaning/disinfection machines..... | 30 |
| | Disposable products..... | 30 |
| 19 | Accessories overview | 32 |
| 20 | Technical specifications | 33 |
| 21 | Signs and symbols | 34 |
| 22 | Warranty and servicing | 36 |
| | Warranty..... | 36 |
| | Servicing/routine check..... | 36 |
| 23 | Disposal | 37 |
| I | Technical documentation / EMC | 366 |
| II | List of accessories | 370 |

Congratulations

By choosing the Vario, you have acquired a high-quality suction pump which is adaptable to your needs. The QuatroFlex suction aggregate is the innovative continuation of the proven Medela piston/cylinder system. As all Medela suction pumps, this pump provides reliable and simple suction. Its simple handling and cleaning as well as the safety features are additional advantages you receive. A comprehensive range of accessories makes the Vario ideally suited to a wide range of medical applications and can be used for continuous operation. Contact us – we will be pleased to advise you.

1 Warnings and safety instructions

WARNINGS

Indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.

CAUTIONS

Indicates a potentially hazardous situation which, if not avoided, could result in minor or moderate injury.

Safety related tip

Indicating useful information about the safe use of the device.

The Vario is approved exclusively for the use as described in these instructions for use. Medela can only guarantee the safe functioning of the system when the Vario is used in combination with the original Medela accessories (collection system, tubings, filters etc. – see annex “II List of accessories”).

**Please read and observe these warning and safety instructions before operation.
These instructions for use must be kept with the device for later reference.**

Please note that these instructions for use are a general guide for the use of the product. Medical matters must be addressed by a physician.

Medela considers herself only responsible for the effect on BASIC SAFETY, reliability and performance of the Vario if it is used in accordance with the instructions for use.

CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician.

Subject to change.

1.1 Warnings and safety instructions for healthcare professionals

WARNINGS

- For use only by medically trained persons who have been adequately trained in suction procedures and in the use of aspirators.
- For patients for whom a breakdown of the device can lead to a critical situation (e.g. patients with acute dyspnoea or severe catarrhal congestion), a replacement device must always be available.
- Gastric suction: When the device is used for gastric suction, the patient's condition and correct functioning of the pump (on the basis of acoustic and visual signs) must be checked at least once every four hours.
- The device must not be used for suctioning explosive, easily flammable or corrosive liquids
- The connecting tubing supplied with the device must never come into direct contact with the suction area. A sterile suction catheter must always be used (risk of infection).
- Before cleaning the device, pull the plug out of the fixed mains socket.
- No modification of this equipment is allowed.
- Do not connect this device to a passive drainage tube.
- Consult the indications for use and consider risk factors and contraindications before using the Vario. Failure to read and follow all instructions in this manual prior to use may result in serious or fatal injury of the patient.
- The Vario pump may shortly shut down with electrostatic discharge (ESD) events at the DC port of 15kV.

CAUTIONS

- Incorrect use can cause pain and injury to the patient.
- Do not use sterile accessories when the sterile packaging is damaged.
- Wireless communications equipment such as wireless home network devices, mobile phones, cordless telephones and their base stations, walki-talkie can affect the Vario and should be kept at least a distance 1ft (30 cm) away from the equipment.
- Portable and mobile RF communications equipment can affect medical devices.
- The patient should be monitored regularly according to the physicians instructions and facility guidelines. Objective indications or signs of a possible infection or complication must be met immediately (e.g. fever, pain, redness, increased warmth, swelling or purulent discharge). Non-observance can lead to considerable danger of the patient. Monitor the Vario frequently for operating status.

Safety related tip

- The Vario suction pump is Magnetic Resonance (MR) unsafe. Do not take the pump into the MR environment.
- For safety tests, it is assumed that the device is serviced and repaired throughout its service life in accordance with the service manual. Because the Vario is a device of safety class II (EN IEC 60601-1), the safety tests are confined to visual inspection of the housing and power cord. This test must be carried out before each use.
- Testing the patient leakage current: Please see service manual for details.

- The protection of the Vario against the effects of the discharge of a cardiac defibrillator is dependent upon the use of appropriate cables.
- Separation from the mains is only assured through the disconnection of the mains plug and the fixed mains socket.
- Before you plug in the device, please check that your local power supply is the same as the voltage given on the specification plate.

Safety instructions

- The Medela Vario is a medical device that requires special safety measures in regard to EMC. It must be installed and put into operation in accordance with the EMC information in annex “I Technical documentation (EMC)”.
- In the case of overflow, inform the internal technical service immediately and perform the tasks in the service manual.
- In each of the following cases, the device must not be used and it must be repaired by Customer Services:
 - if the power cord or the plug are damaged
 - if the device is not functioning perfectly
 - if the device is damaged
 - if the device shows clear safety defects.
- Keep the power supply cord away from hot surfaces.
- The mains plug and the on/off switch must not come into contact with moisture.
- Never pull the mains plug out of the fixed mains socket by pulling on the power supply cord!
- Never leave the device unattended when it is switched on.
- The pump must stand upright during use.
- Never place the power supply cord around your neck.
- Never use the device at high room temperatures, while bathing or showering, if you are very tired or in an environment where there is a risk of explosion.
- Never place the device in water or other liquids.
- When using single use, sterile products, please note that they are not intended to be reprocessed. Reprocessing could cause loss of mechanical, chemical and/or biological characteristics. Reuse could cause cross contamination.
- Contact your local Medela customer service representative for assistance with product operations.
- Do not place a 2l suction jar in the lateral adaptor.
- The overflow protection filter protects the pump against overflow. Caution: Suctioning is interrupted when the filter is wet/moist or clogged. Test the filter before each use and during use periodically (see chapter 17).

These instructions for use must be kept for later reference.

1.2 Warnings and safety instructions for patients and caregivers

Read the entire instruction manual before trying to operate the device.

If you have any questions about your Vario suction pump, please contact your healthcare provider/physician.

WARNINGS

- The individual responsible for operating this device on a daily basis (patient, caregiver) must be able to read, understand and follow directions provided in these instructions for use, and by the healthcare professional. The individual must be able to hear and respond to auditory signal.
- Only use the Vario suction pump on the person for whom it was ordered and only for its intended use. Do not perform the therapy without your healthcare provider's supervision.
- If you are dependent on the device for airway suctioning and a breakdown can lead to a critical situation, you must have a replacement device at hand.
- The connecting tubing delivered with the device must always be connected to a sterile catheter as prescribed by your physician.
- Do not modify the equipment.
- Vario is EMC-tested in conformity with the requirements of IEC 60601-1-2:2007 and IEC 60601-1-2:2014 4th Edition according to clause 7 and 8.9. Vario is suitable for use in home and clinical environments.
- Wireless communications equipment such as wireless home network devices, mobile phones, cordless telephones and their base stations, walki-talkie can affect the Vario and should be kept at least a distance 1ft (30 cm) away from the equipment.

General safety at home

- Vario should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, Vario should be observed to verify normal operation in the configuration in which it will be used.
- Keep the power supply cord away from hot surfaces.
- Keep the mains plug and the on / off switch away from moisture.
- Never pull the mains plug out of the fixed mains socket by pulling on the power supply cord!
- Never leave the device unattended when it is switched on.
- The pump must stand upright during use.
- Never use the device at high room temperatures, while bathing or showering or if you are very tired.
- Do not use extension cords with Vario pumps. Keep the power supply cord where your will not fall or trip over it.
- Never place the power supply cord around your neck.
- Keep the Vario pump clean and dry.

- Never place the pump in water or liquids.
- If the pump gets wet, rub with dry towel. Do not dry in a microwave.
- Keep the pump away from children and pets.
- Keep the pump protected from direct sunlight.
- Prevent the pump from falling to the floor during sleep.
- Do not touch the pump when it has fallen into water. Unplug device immediately from electrical outlet.
- The Vario pump may shortly shut down with electrostatic discharge (ESD) events at the DC port of 15kV.

When you leave the house

- Make sure the LED indicator is not blinking when you leave the house.
- For your convenience put the pump into the optional carrying bag and store the tubing safely in the tubing bag.
- Don't forget to take a spare jar/liner/tubing and the power supply with you.

These instructions for use must be kept for later reference.

2 Power supply and battery operation

The Vario is a mains-powered suction pump. Before you plug in the device, please check that your local power supply is the same as the voltage given on the specification plate. If the Vario AC/DC versions are running continuously, a maximum of 30 minutes of battery time is expected before it is necessary to recharge. If the Vario AC/DC version is not in use, the battery must be calibrated every 60–90 days.

Battery low (only AC/DC versions)

During battery operation a slow blinking of the LED indicator and an acoustic signal (beep - beep - beep) indicates, that the rechargeable battery is soon discharged. The function of the Vario is maintained for 10 minutes (20 minutes Vario low vacuum versions), but the battery must be recharged as soon as possible.

Battery empty (only AC/DC versions)

As soon as the rechargeable battery is empty, the pump comes to a stop. The proper function of the pump is no longer maintained. The rechargeable battery must be recharged.

3 Description

Introduction

The Vario is a high-quality suction pump. It is powered by the well-proven QuatroFlex system and guarantees maximum suction performance for many suctioning needs. It ideally combines easy handling and cleaning with safety features to ensure optimal operation. You can choose from a comprehensive range of accessories from Medela to configure the pump to many medical applications. The device should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the device should be observed to verify normal operation in the configuration in which it will be used.

Intended use

The Vario 8/18/ci suction pumps are indicated for aspiration and removal of surgical fluids, tissue (including bone), gases, bodily fluids (including vomit) or infectious materials from a patient's airway or respiratory support system, either during surgery or at the patient's bedside. Generally the Vario 8/18/ci is intended to be used for a variety of suctioning procedures including nasopharyngeal, tracheal, surgical, gastrointestinal suctioning either in «constant» or «intermittent» mode. Please make sure you have a correct Vario version for the required application.

Contraindications

None known.

Intended user

The Vario should only be operated by properly trained staff or users (e.g. caregivers, patients). These persons must not be hard of hearing or deaf and must have adequate visual faculty. The training has to be refreshed at least once a year.

Intended Patient Population

The Vario is intended to be used on patients only exhibiting conditions as described in the indications for use.

Important note

Compliance with proper surgical procedures and techniques is the responsibility of the physician. Each physician must evaluate the appropriateness of the treatment based on his own knowledge and experience.

Expected service life

The expected service life of the device is seven years.

4 Overview

Definition of vacuum

By the application of medical aspiration devices, vacuum is normally given as the difference (in absolute figures) between absolute pressure and atmospheric pressure or as negative values in Kilopascal (kPa). In this document, the indication of -10 kPa for example always refers to a pressure range in kPa below atmospheric ambient pressure (according to EN ISO 10079-1).

Versions of the pump

Vario 8 AC

Vario 8 AC/DC

Vario 8 c/i AC/DC

Vario 18 AC

Vario 18 AC/DC


Vario 18 c/i AC/DC


Front of the pump

Bottom of device

Operating elements and battery status

not lit Battery is fully charged

lit up – Battery is being charged
– Pump is connected to the mains


blinks slowly Battery is soon empty

blinks rapidly Battery is defective


Acoustic indicator

During battery operation an acoustic signal (beep - beep - beep) indicates that the rechargeable battery is soon discharged.


On/off switch AC versions

On/off switch AC/DC versions

On/off switch c/i versions

5 Installation

Check initial delivery

Check the delivery package of the Vario for completeness and general condition.

Vario portable version

| | | |
|------------|--------------------|----------|
| REF | Vario 8 AC | 026.XXXX |
| | Vario 8 AC/DC | 026.XXXX |
| | Vario 8 c/i AC/DC | 026.XXXX |
| | Vario 18 AC | 026.XXXX |
| | Vario 18 AC/DC | 026.XXXX |
| | Vario 18 c/i AC/DC | 026.XXXX |

Silicone Tubing \varnothing 7 x 12 mm, 25 cm with 2 coupling pieces

REF 077.0931

5x Disposable overflow protection filter with Medela connections

Instructions for use

REF 177.2604

Accessories

Please ask your healthcare provider/physician in case of missing parts or for additional accessories.

6 Preparation for use

WARNINGS

For use only by persons who have been adequately trained in suction procedures and in the use of aspirators.

CAUTIONS


- The Vario must remain in an upright position during use.
- Sterile accessories must be checked on the integrity of the packaging before use.
- Non sterile and reusable accessories must be cleaned, disinfected and/or sterilized according to chapter 18

6.1 Checks before use

- Check the Vario system before use for damage of the power cord or plug, obvious device damage or safety defects and proper functioning of the device.
- Check for completeness and general condition of the Vario delivery package.
- Only for AC/DC versions: make sure that the rechargeable battery is calibrated according to chapter 16.
- Check all accessories prior to use:
 - suction jars, lids and liners for cracks, brittle and flawed spots. Replace if necessary.
 - Tubing for cracks, brittle areas and that connectors are firmly attached. Replace if necessary.
 - As an additional safety test, evacuate the system (including jars) to maximum vacuum several times before actual use.

6.2 Assembly of the basic configuration

1.

- 1.1 Attach the overflow protection filter to the Vario with the arrow pointing in the flow direction.
For home use: The filter must be exchanged at least on a weekly basis.

2. Attach all necessary accessories according to your needs. See annex "II List of accessories".
For patients and caregivers: Align with your healthcare provider/physician on the accessories you need.

WARNING

The filter protects the pump against overflow.

When the filter is wet/moist/clogged you will not be able to suction.


Test the filter before each use (see chapter 17 «filter test»).

6.3 Assembly of the REUSABLE collection system


positions for jars:

1. If you use a DISPOSABLE collection system proceed with step 6.4.

2.1 Attach mechanical overflow protection to lid.
2.2 Pull gently downwards to make sure it is open/deactivated.

3.1 Attach the lid to the jar.
3.2 Lock with the two lid clamps.

4.1 Attach the jar to the Vario.
4.2 Connect tubing to the filter.
4.3 Connect tubing to the lid of the jar (marked with «vacuum»).

5.1 Connect patient tubing to the lid of the jar (marked with «patient»).

CAUTION

Do not use tubing or other sterile accessories if the sterile packaging is damaged.
Do not reuse single use or sterile products.

6.4 Assembly of the DISPOSABLE collection system

positions for jars:


1. If you use a REUSABLE collection system proceed with step 6.3.


2.1 Attach the Vario jar adapter to the Vario.

3.1 Choose the 1.5l size (yellow code) for the Vario.
3.2 Prepare and insert corresponding liner size (yellow) according to instructions delivered with liners. The disposable liner has an integrated overflow protection filter. No action is necessary to activate it.

4.1 Attach the jar to the Vario.
4.2 Connect tubing to the filter.
4.3 Connect tubing to the jar. The tubing connector is on top of the jar.

5.1 Connect patient tubing to the lid of the liner using either the coloured angle piece or the patient port directly (depending on tubing size).

CAUTION

Do not use tubing or other sterile accessories if the sterile packaging is damaged.
Do not reuse single use or sterile products.

7 Operating instructions

CAUTION

The Vario is to be set up in such a way, that a separation from the mains supply can be easily managed.

7.1 Connect Vario to mains power

1. Check the pump before use following the instruction in chapter 6.1.

2. **All versions: if the pump is connected to a fixed mains socket:**

Plug in the mains plug of the power cord to a fixed mains socket.

or

AC/DC versions: if the pump is operated with a 12VDC power source:


Connect the optional 12VDC connection cable to the socket on the side of the pump and connect it to the 12VDC power source.

or

AC/DC versions: if the pump is operated with the rechargeable battery:

No connection to a power source is needed. Make sure that the battery was charged for five hours.

7.2 Check maximal vacuum for functional check

1.1 Switch on the Vario. The pump will start running. Caution: **c/i versions** must be switched on in «constant» mode.

2.1 Press and turn the vacuum regulator to the right to set maximum vacuum.

3.1 Seal the end of the patient tubing with your thumb.
3.2 Compare the maximum vacuum according to the specification. See chapter 10 if the maximum vacuum is not reached.

Specifications:

| Altitude above sea level: | Vario 8 / Vario 8 c/i Max. Vacuum: low vacuum | Vario 18 c/i Max. Vacuum: medium vacuum | Vario 18 Max. Vacuum: high vacuum |
|---------------------------|---|---|---|
| + 2000 m | - 7.0 kPa | - 43 kPa | - 59 kPa |
| + 1000 m | - 7.9 kPa | - 48 kPa | - 66 kPa |
| + 500 m | - 8.5 kPa | - 51 kPa | - 70 kPa |
| 0 m | - 9.0 kPa | - 55 kPa | - 75 kPa |


(Tolerance +/- 10%)

7.3 Changing vacuum level


- 1.1 Clamp patient tubing
- 1.2 Push and turn vacuum regulator to select the correct vacuum according to the particular application.
- 1.3 Check vacuum gauge for setting.

7.4 Placing out of operation after use

- 1.1 Switch the on/off switch of the Vario to position «0» (c/i - versions: «Off»).
- or
2. **If the pump is connected to a fixed mains socket:** Disconnect the mains plug from the fixed mains socket.
 - or
 - If the pump is connected to a 12VDC power source:** Disconnect the 12VDC plug from the 12VDC power source.
3. Clean and disinfect the Vario. See chapter 18.

8 Changing reusable jar

- 1.1 Switch the on/off switch of the Vario to position «0» (c/i - versions: «Off»).
- 1.2 Remove patient and vacuum tubings from the lid.

- 2.1 Remove full jar from the Vario.

- 3.1 Prepare new jar and unpack it.
- 3.2 Secure new jar to the Vario
- 3.3 Reconnect vacuum tubing and new patient tubing firmly.

- 4. Empty jar and dispose of patient tubing in accordance with local guidelines and inhouse guidelines.
- 5. Switch on the Vario. The pump will start running.
- 6. Set vacuum according to the particular application. See chapter 7.3.


9 Changing disposable liner


WARNING

Disposable liners are not intended to be reprocessed. Reprocessing could cause loss of mechanical, chemical and/or biological characteristics.

Single use only. After use, dispose of according to local regulations. Do not clean or reprocess. Reuse can lead to mechanical failure of the liner can lead to exposure to biological hazards of patient and staff.

- 1.1 Remove patient tubing with colored angle piece from lid.
- 1.2 Plug patient port.


- 2.1 Switch the on/off switch of the Vario to position «0» (c/i - versions: «Off»). Reduce vacuum to minimum by pressing and turning vacuum regulator to the left.

- 3.1 Remove full liner from the jar.

4. Dispose of the liner and patient tubing in accordance with local guidelines and inhouse guidelines.
5. Prepare new liner, unpack it and insert it in the jar. Press down lid.

- 6.1 Connect new patient tubing firmly.

7. Switch on the Vario. The pump will start running.
8. Set vacuum according to the particular application. See chapter 7.3.

10 Troubleshooting

Motor not running

Check if:

- the Vario is switched on.
- the mains voltage is correct and if the mains plug is inserted correctly into the fixed mains socket.
- the internal battery is charged (only AC/DC version under battery operation)
- the fuse on the bottom of the Vario is not defective. For replacing the defective fuse see chapter 11.

Insufficient vacuum

Check if:

- the vacuum regulator is set correctly.
- the tubings are not defective or broken. If necessary, replace.
- all plug-in connections are tight.
- the floater on the reusable lid moves freely up / down.
- the overflow protection is deactivated/open. If the overflow protection is activated, deactivate as shown under 6.3a/2.2.
- the suction jar/lid have no cracks, brittle areas, discolouration. Replace if necessary.
- the disposable system has no cracks, brittle areas, discolouration. Replace if necessary.
- filter is not clogged. To test if the filter is clogged see chapter 17.
- the QuatroFlex is defective. To replace the QuatroFlex see chapter 13.

No LED lit (AC/DC version only)

The rechargeable battery is fully charged.

LED lit up (AC/DC version only)

- The rechargeable battery is being charged.
- The pump is connected to the mains.

LED blinks slowly (AC/DC version only)

The rechargeable battery is soon empty.

LED blinks rapidly (AC/DC version only)

The rechargeable battery is defective. For replacing the defective battery see chapter 14.

CAUTIONS

Do not try to fix or repair the Vario pump.

Contact your healthcare provider/biomed immediately if:

- if an overflow into the device has occurred
- if the power cord or the plug are damaged
- if the device is not functioning perfectly
- if the device is damaged
- if the device shows clear safety defects.

If the fault cannot be rectified, please contact your healthcare provider/biomed.

11 Replacing defective fuse

Only for Healthcare professionals

WARNINGS

Before replacing the fuse, pull the mains plug from Vario out of the fixed mains socket.

1. Disconnect the Vario from the mains.

- 2.

- 2.1 Open (turn left) the fuse holder on the bottom of the Vario.

3. Replace the defective fuse(s). Make sure that the technical specification are corresponding.

AC-versions: **T 0.8 AL, 5 x 20, 230–240 V**

T 1.25 AL, 6.3 x 20, 120 V

AC/DC-versions: **T 1.0 AL, 5 x 20, 100–240 V**

- 4.

- 4.1 Close (push back and turn right) the fuse holder.

5. Reconnect the Vario to the mains again.

6. Switch on the Vario again.

12 Setting interval times (c/i versions only)


Only for Healthcare professionals

WARNINGS

Before setting the interval times, switch off the Vario and pull the mains plug from Vario out of the fixed mains socket.


1.

1.1 Unscrew the 4 screws on the back of the Vario.

1.2 Open the Vario by removing the back cover.

2.

2.1 Determine the location of the red jumper.

2.2 Place the red jumper according to the desired interval:

| On-time | Off-time |
|----------------|-----------------|
| 8 sec. | 4 sec. |
| 16 sec. | 8 sec. |
| 32 sec. | 16 sec. |

3. Reattach the back cover to the Vario and tighten the 4 screws.

4. Plug the Vario to a fixed mains socket and test the correct functioning of the Vario (see chapter 7).





13 Replacing the QuatroFlex

Only for Healthcare professionals

WARNINGS

Before replacing the QuatroFlex, switch off the Vario and pull the mains plug from Vario out of the fixed mains socket.

1. 
 - 1.1 Unscrew the 4 screws on the back of the Vario.
 - 1.2 Open the Vario by removing the back cover.
2. 
 - 2.1 Remove the 3 marked tubings from the QuatroFlex.
3. 
 - 3.1 **Turn** the QuatroFlex 45° to the left.
 - 3.2 **Pull out** the QuatroFlex and remove it from the Vario.
4. Clean inside of the Vario.
5. 
 - 5.1 **Insert** new QuatroFlex.
 - 5.2 **Turn** it 45° to the right.
 - 5.3 **Reattach** the 3 tubings to the QuatroFlex.
6. Reattach the back cover to the Vario and tighten the 4 screws.
7. Plug the Vario to a fixed mains socket and test the correct functioning of the Vario (see chapter 7).

14 Replacing the rechargeable battery (AC/DC versions only)


Only for Healthcare professionals

WARNINGS

Before replacing the rechargeable battery, switch off the Vario and pull the mains plug from Vario out of the fixed mains socket.


1.

1.1 Unscrew the 4 screws on the back of the Vario.


1.2 Open the Vario by removing the back cover.

2.

2.1 Detach the battery cable by pulling the plug.

3.

3.1 Remove the old rechargeable battery.

3.2 Insert the new rechargeable battery.

4.

4.1 Reattach the battery cable.

5. Reattach the back cover to the Vario and tighten the 4 screws.

6. Plug the Vario to a fixed mains socket and calibrate the rechargeable battery (see chapter 16).

15 Battery data and test (AC/DC versions only)

CAUTIONS

Only charge the rechargeable battery when the acoustic signal sounds and the charging indicator blinks slowly. The rechargeable battery is charged during mains operation.

The Vario 18 is a mains- and battery powered suction pump. One full battery charge (Vario 18 AC/DC version) would last you a total of 30 minutes continuous running time. Then the battery must be recharged.

If you are not using your Vario 18 AC/DC version or storing for a longer period of time (> 2 months), recharge the battery every 60–90 days.

Battery low

Only recharge the battery when the LED indicator blinks slowly and you hear the acoustic signal (beep - beep - beep). You have 10 minutes of suctioning time left.

Battery empty

As soon as the rechargeable battery is empty, the pump comes to a stop. The rechargeable battery must be recharged.

Duration of pumping

If the battery is fully charged, the operating time is about:


- low vacuum versions: approx. 120minutes
- high and medium vacuum versions: approx. 30minutes.

Charging time (charged up to 100%)


- low vacuum versions: approx. 5hours
- high and medium vacuum versions: approx. 5hours.

Testing the rechargeable battery

1. Make sure that the rechargeable battery is fully charged.
2. Pull out the mains plug from the fixed mains socket.

- 3.1 Press vacuum regulator and turn it to the right to select maximum vacuum.

- 4.1 Switch on the Vario. The pump will start running. Caution: **c/i versions** must be switched on in «constant» mode.

5. Measure the time until the pump comes to a stop. The battery is ok, when the following values are reached:
 - low vacuum versions: > 120 min.
 - high and medium vacuum versions: > 30 min.

If the values are not reached, repeat test or calibrate the battery (see next page).

16 Battery calibration (AC/DC versions only)

CAUTIONS


Calibrate rechargeable batteries before the first use, or after a storage period >2 months. Replace the battery if the calibration process fails.

Calibrating NiMH rechargeable battery

1. Make sure that the battery is fully charged.

2. Pull out the mains plug from the fixed mains socket.

3.

3.1 Switch on the Vario. The pump will start running. Caution: **c/i versions** must be switched on in «constant» mode.

4. Let the pump run until the rechargeable battery is empty and the pump comes to a stop.

5. Plug the mains cable into a fixed socket. The rechargeable battery is being charged.

6. Wait until the charging indicator goes out.


7. Repeat step 1 to 6. The Vario AC/DC is then ready to use.


Storage of the rechargeable battery


In order to keep the self-discharge at a minimum, store the Vario and the replacement rechargeable batteries at temperatures below 25°C (77°F). Repeat the calibrating process every 60–90 days.

17 Filter test

1.  1.1 Remove tubing from filter.

2.  2.1 Select maximum vacuum.

3.  3.1 Switch on the Vario. The pump will start running. Caution: **c/i versions** must be switched on in «constant» mode.

4.  4.1 Read the vacuum.
4.2 Replace the filter if the vacuum exceeds the values given in the following table:

| | |
|--|---------|
| Vario 8 / Vario 8 c/i low vacuum | - 3kPa |
| Vario 18 high vacuum | - 10kPa |
| Vario 18c/i medium vacuum | - 10kPa |

18 Cleaning guidelines

WARNINGS

After each use, the parts that had been in contact with the aspirated secretions are to be cleaned, disinfected, sterilised or disposed of **according to the table on the next page**.

Before cleaning the device, pull the mains plug out of the fixed mains socket.

General notes

- These are general recommendations only that may be adjusted individually, based on the hospital's specific directives and cleaning practices and policies.
- Detailed information according to Medela cleaning instructions (product code 200.2391).
- Wear protective gloves for cleaning/disinfection.
- Dispose of fluids such as blood and secretions and the parts contaminated with them in accordance with internal hospital guidelines.

Medela recommended surface cleaning agents for pump housing

- CaviWipes, Metrex Research, www.metrex.com
- Mikrozid AF Wipes, Schülke & Mayr, www.schuelkemayr.com

Water

Use only the purest quality of water for cleaning. Water hardness is a serious consideration since deposits left on medical products may not be properly decontaminated. Use deionised water in order to reduce this problem. The final rinse water is bacterial free and contains no endotoxins.

Cleaning/disinfection machines

Can be used to desinfect parts from the table on the next page. A hot water rinse (maximum temperature 100 °C) may provide a medium-to high level of disinfection. Every section of the constituent parts must be accessible in order to ensure efficient cleaning. We recommend using a cleaning/disinfecting machine that has been approved by the Robert Koch- Institute and complies with ISO 15883. Recommended temperature for noncritical medical devices (i.e. those that only come into contact with uninjured skin) is 90 °C for 1 minute. The time is increased to 5 minutes for all medical devices that are considered to be critical.

Disposable products

These are single use products not intended to be reprocessed. Reprocessing could cause loss of mechanical, chemical and/or biological characteristics. Reuse could cause cross contamination.

| PSU* jars, lids | PC** jars | Silicone tubing | Pump housing Mains cable | Plastic accessories | Legend: * Polysulfone ** Polycarbonate |
|-----------------|-----------|-----------------|-----------------------------|---------------------|--|
| x | x | x | x | x | 1. Disassembly Separate all individual parts before cleaning, disinfecting and sterilising. |
| x | x | x | | | 2. Cleaning Clean components in hot water (60–70 °C) containing a detergent with a pH range between 6.0 and 8.0 only, in order to avoid damaging the instruments and containment devices. Enzymatic detergents help with the removal of organic matter, such as blood. Detergents should be used according to their manufacturer's recommended concentration levels. Some alkaline detergents have been formulated to be safe for reprocessing medical devices. The manufacturers of these detergents should provide information about specific materials that may be damaged by their detergents. Soak all parts thoroughly with warm, soapy water or in enzymatic detergent for 1–5 minutes. 1. Remove visible dirt with a cleaning tool – general purpose cleaning brushes, such as pipe cleaners or non-abrasive lint cloths. Brushes and pipe cleaners should fit snugly but still be able to be moved around easily within the area to be cleaned. Rinse thoroughly in clear water. Allow to dry 2. Check the instruments for visible dirt and repeat these steps if necessary. |
| | | | x | x | 3. Cleaning Wipe with detergent mentioned above. Use clean cloth to dry. |
| x | x | x | | x | 4. Disinfect Soak at room temperature for 30 minutes in a disinfection/ cleaning solution (A). After the reaction time, use water to rinse the solution residue from the individual parts and rinse the tubing. Rinse for at least 2 minutes with the cleaning solution (B). (A) fluid aldehyde-free disinfection solution with cleaning effect for use as a soaking bath, bactericide, fungicide, limited virucide, with good material compatibility for stainless steel, non-ferrous material and plastics including silicone, slightly alkaline. e.g. 2 % neodisher® Septo MEDsolution (B) Fluid, pH-neutral, enzymatic cleaning solution for the treatment of instruments by machine or manually with very good material compatibility for stainless steel, non-ferrous material and plastics including silicone. or use cleaning/ disinfecting machine. or use your dishwasher at home. |
| x | | x | | | 5. Sterilise Remove any chemical disinfectant residue before autoclaving. Do not stack during autoclaving. In the autoclave at 134 °C for 10 minutes. The sterilisation device must comply to ISO 17665:2003. The following procedures are acceptable: In the autoclave with saturated steam at 132 °C, with triple pre-evacuation, sterilisation time of 6 minutes, or at 135–137 °C in an porous load vacuum autoclave for 3–3.5 minutes |
| x | x | x | | x | 6. Storing components After disinfection/ sterilising, store the components in sterile foil until required for use. |

19 Accessories overview

WARNINGS

Vario was verified in combination with the accessories listed in annex “II List of accessories”. For a correct and safe operation use Vario with these accessories only. Further information is supplied on the instruction sheet of the individual accessory.

Safety related tip

If the pump is used together with non-Medela applied parts, they must:

- have CE mark and (if necessary) local registration
- be able to be attached to Medela accessories safely without impacting the performance of the pump.

Note: when combining Medela parts and non-Medela parts, you take over the responsibility of the entire system.

20 Technical specifications

CAUTION

Do not operate the Vario 18 Pump above 9,842 feet above sea level (3000 m above sea level).

Vario 8/Vario 8 c/i:

low vacuum, -9 kPa/-68 mmHg (Tolerance: ±10%)

Vario 18 c/i:

medium vacuum, -55 kPa/-413 mmHg (Tolerance: ±10%)

Vario 18:

high vacuum, -75 kPa/-563 mmHg (Tolerance: -10%)

Measured at 0m, atmospheric pressure: 1013.25 hPa
Please note: vacuum levels may vary depending on location (meters above sea level, atmospheric pressure and temperature).

Vario 8/Vario 8 c/i: 8 l/min. (±10%)

Vario 18/Vario 18 c/i: 18 l/min. (±10%)

AC 3.5 kg
AC/DC 4.2 kg
Without jar

230-240V, 50 Hz, 90 VA
230-240V, 60 Hz, 90 VA
120V, 60 Hz, 70 VA

100-240V
50/60 Hz
80 VA

Car adapter
12 VDC / 3A
8A (Fuse)


ISO 13485
CE (93/42/EEC), IIa

0123

h x w x d
380 x 170 x 285 mm


Transport / Storage Conditions

The Vario pump and accessories must remain in the packaging for storage and stored at a temperature range of -25 °C and +70 °C (-13 °F and +122 °F). Rechargeable battery 0 °C and +25 °C (+32 °F and +77 °F).

Transport / Storage

Rechargeable battery

Operating Conditions

The Vario Pump and accessories must be operated within a temperature range of +5 °C and +40 °C (+41 °F and +104 °F). Do not operate the products in extreme cold or heat. Avoid high humidity (15-93 % is permissible).

IP21

The Vario Pump is protected against ingress of dripping water.

21 Signs and symbols

This symbol indicates the compliance with the essential requirements of the Council Directive 93/42/EEC of 14 June 1993 concerning medical devices.

This symbol indicates the class of the pump.

This symbol indicates a class II device.

IP21

This Symbol indicates the protection against ingress of solid foreign objects and against harmful effects due to the ingress of water.

This symbol indicates that the device should not be used after the date shown.

This symbol indicates a single use device. Do not reuse the device.

REF

This symbol indicates manufacturer's catalogue number.

SN

This symbol indicates manufacturer's serial number.

LOT

This symbol indicates manufacturer's batch code.

STERILE EO

This symbol indicates the device is sterilized using ethylene oxide.

This symbol indicates the temperature limitation for operation, transport and storage.

This symbol indicates the humidity limitation for operation, transport and storage.

This symbol indicates the compliance with additional USA and Canada safety requirements for medical electrical equipment.

This symbol indicates a type CF applied part.

This symbol indicates MR UNSAFE.

This symbol indicates that interferences may occur in the vicinity of equipment marked with this symbol.

This symbol indicates the manufacturer.

This symbol indicates the date of manufacture.

This symbol indicates the legal specifications of the pump.

This symbol indicates to follow instructions for use.

This symbol indicates to consult instructions for use.

This symbol indicates a CAUTION or WARNING associated with the device.

This symbol indicates safety related tip.

This symbol indicates the atmospheric pressure limitation for operation, transport and storage.

This symbol indicates do not use the device if package is damaged.

This symbol indicates the number of items n that the content is sufficient for.

This symbol indicates do not dispose the device together with unsorted municipal waste (for EU only).

This symbol indicates that the material is part of a recovery/recycling process.

This symbol indicates a carton package.

This symbol indicates to keep the device away from sunlight.

This symbol indicates to handle the fragile device with care.

This symbol indicates to keep the device dry.

This symbol indicates the maximum vacuum level of the pump.

This symbol indicates the flow levels of the pump.

This symbol indicates the electrical specifications of the pump.

This symbol indicates alternating current.

This symbol indicates direct current.

This symbol indicates the weight of the pump.

This symbol indicates the dimensions (h x w x d) of the pump.

This Symbol indicates a Prescription Device. CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician. (for US only).

This symbol indicates number of items.

This symbol indicates the fuse

This symbol indicates the item is a medical device.

This symbol indicates the authorized representative in the European Community/ European Union

This symbol indicates the quantity (X) of individual devices in pack.

This symbol indicates the specifications to connect to car adapter.

22 Warranty and servicing

Warranty

Medela AG warrants the device will be free from defects in materials and workmanship for a period of 2 years (6 month for the rechargeable battery and the QuatroFlex) from the date of delivery. Faulty material will be replaced free of charge during this period if not resulting from abuse or misapplication. This will not apply to parts subject to wear and tear in use. To ensure compliance with this warranty as well as optimum service from Medela products, we recommend the exclusive use of Medela accessories with our appliances.

In no event shall Medela AG be liable for claims which exceed the scope of warranty described including liability for consequential damages, etc. The right to the replacement of faulty parts will not be recognized by Medela if any work has been carried out on the pump by unauthorized persons. This warranty is subject to the appliance being returned to a Medela service centre.

Servicing/routine check

Routine checks and service work are only to be carried out by positions authorised by Medela. Medela recommends a routine check every year (according to product specific Service Manual). The Medela Service Manual will be made available for authorized positions.

23 Disposal

The Vario comprises metals and plastics and should be disposed of in accordance with the European directives 2011/65/EU and 2012/19/EU. Additional, local guidelines must also be observed. With the AC/DC versions, the electronic components and the rechargeable battery must be disposed of separately, in accordance with the local regulations.

Please take care that you dispose the Vario and its accessories in accordance with your local disposal guidelines.

For patients and caregivers : Ask your healthcare provider how to dispose of used accessories (e.g. tubings, filters, liners).

User information for the disposal of electrical and electronic equipment

This symbol means that the electrical and electronic equipment must not be disposed as normal household refuse. A correct disposal of this device protects and prevents possible damage to the environment or human health. For more information about the disposal contact the manufacturer, your local caregiver or healthcare provider. This symbol is only valid in the European Union. Please respect the relevant state laws and rules in your country for the disposal of electrical and electronic equipment.

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I Technical documentation

WARNINGS

Do not use other accessories than those specified or sold by the manufacturer as replacement parts for internal components as it may result in increased emissions or decreased immunity of the Vario pump.

HF (high-frequency) surgical equipment, radio networks or the like can influence the operation of the device and may not be operated in combination with the system.

EMC

Vario is EMC-tested in conformity with the requirements of IEC 60601-1-2:2007 and IEC 60601-1-2:2014 4th Edition according to clause 7 and 8.9. Vario is suitable for use in home and clinical environments. Vario is a medical device that requires special safety precautions and must be installed and placed in operation in accordance with the attached EMC information.

Electromagnetic Compatibility (EMC, IEC 60601-1-2:2007, 3rd Edition and IEC60601-1-2 4th Edition, Table 1)

Electromagnetic Emissions

Vario is intended for use in the electromagnetic environment specified below. The customer or the user of Vario should assure that they are used in such environment.

| Emission Tests | Compliance | Electromagnetic environment – guidance |
|---|----------------|--|
| RF emissions CISPR 11 | Group 1 | Vario uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference with nearby electronic equipment. |
| RF emissions CISPR 11 | Class B | Vario is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| Harmonic emissions IEC 61000-3-2 | Class A | |
| Voltage fluctuations/ flicker emissions IEC 61000-3-3 | Not applicable | |

WARNINGS

Vario should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, Vario should be observed to verify normal operation in the configuration in which it will be used.

Electromagnetic Immunity

Vario is intended for use in the electromagnetic environment specified below. The customer or the user of Vario should assure that they are used in such environment.

| Immunity Tests | IEC 60601 test level | Compliance level | Electromagnetic environment – guidance |
|--|---|---|--|
| Electrostatic Discharge (ESD) IEC 61000-4-2 | ± 6 kV contact ± 8 kV air | ± 6 kV contact ± 8 kV air | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%. |
| Electrical fast transient/burst IEC 61000-4-4 | ± 2 kV for power supply lines ± 1 kV for input / output lines | ± 2 kV for power supply lines ± 1 kV for input/ output lines | Mains power quality should be that of a typical commercial or hospital environment. |
| Surge IEC 61000-4-5 | ± 1 kV line(s) to line(s) ± 2 kV line(s) to earth | ± 1 kV line(s) to line(s) ± 2 kV line(s) to earth not applicable | Mains power quality should be that of a typical commercial or hospital environment. |
| Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11 | <5% U_T (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec | <5% U_T (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec | Mains power quality should be that of a typical commercial or hospital environment. If the user of Vario requires continued operation during power mains interruptions, it is recommended that Vario be powered from an uninterruptible power supply or a battery. |
| Power frequency (50/60 Hz) magnetic field IEC 61000-4-8 | 3 A/m | 3 A/m | Power frequency magnetic fields should be at levels of a typical commercial or hospital environment. |

Note: U_T is the a.c. mains voltage prior to application of test level.

Electromagnetic Immunity

Vario is intended for use in the electromagnetic environment specified below. The customer or the user of Vario should assure that Vario is used in such environment.

| Immunity Tests | IEC 60601 test level | Compliance level | Electromagnetic environment – guidance |
|-------------------------------|--------------------------------|------------------|--|
| Conducted RF IEC 61000-4-6 | 3 Vrms 150 kHz to 80 MHz | 3 Vrms | <p>Portable and mobile RF communications equipment should be used no closer to any part of Vario, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = 1.2 \sqrt{P}$ $d = 1.2 \sqrt{P}$ <p>80 MHz to 800 MHz</p> $d = 2.3 \sqrt{P}$ <p>800 MHz to 2.5 GHz</p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m)</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range.^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p>  |
| Radiated RF IEC 61000-4-3 | 3 V/m 80 MHz to 2.5 GHz | 3 V/m | |

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed RF transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which Vario are used exceeds the applicable RF compliance level above, Vario should be observed to verify normal operation. If abnormal operation is observed, additional measures may be necessary, such as reorienting or relocating Vario.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distance between portable and mobile RF communications equipment and Vario

Vario is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of Vario can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and Vario as recommended below, according to the maximum output power of the communications equipment.

| Rated maximum output power of transmitter | Separation distance according to frequency of transmitter | | |
|---|---|---|--|
| | 150 kHz to 80 MHz $d = 1.2 \sqrt{P}$ | 80 MHz to 800 MHz $d = 1.2 \sqrt{P}$ | 800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$ |
| W | | | |
| 0.01 | 0.12 | 0.12 | 0.23 |
| 0.1 | 0.38 | 0.38 | 0.73 |
| 1 | 1.2 | 1.2 | 2.3 |
| 10 | 3.8 | 3.8 | 7.3 |
| 100 | 12 | 12 | 23 |

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

II List of accessories

Patient tubing (from collection system to patient applied part)

| | | | |
|---|---|---|---|
| Disposable, 7x10 mm, sterile | | Silicone, 7x12 mm | |
|  |  |  |  |
| | 077.0170: 150 cm | 077.0053: 32 cm | 077.0055: 120 cm |
| | | 077.0055: 120 cm | 077.0056: 150 cm |
| | | 077.0056: 150 cm | 077.0051: 200 cm |
| | | 077.0051: 200 cm | 077.0164: 300 cm |
| | | 077.0164: 300 cm | 077.0052: 2500 cm |
| | | 077.0052: 2500 cm | 077.0961: 180 cm, silicone (green) |
| Disposable, 6x8 mm, sterile | | Disposable, 12x18 mm | |
|  |  |  |  |
| | 077.0184: 150 cm | | 077.0963: 180 cm |
|  |  | Silicone, 5x10 mm | |
| | 077.0965: 150 cm |  |  |
| Disposable, 6x8 mm | | | 077.0970: 2500 cm |
|  |  |  | 077.1017 |
| | 077.0951: 180 cm, sterile |  | 077.1020 |
|  |  | | |
| | 077.0952: 180 cm | | |
|  |  | | |
| | 077.0191: 300 cm, sterile | | |

*

Liners

077.0083/84
101033142/
101032951

Reusable lids

077.0440 077.0420/30

PC jars (max. 4x2.5l)

077.0082 077.0105

Reusable jars (max. 4x5l)

077.0125: 0.25l
077.0155: 0.5l
077.0110: 1l


Reusable sets

077.0715: 1l

Vacuum tubing (from suction pump to collection system)

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Medela UK Ltd.
Medela LLC